

Product datasheet for SC335918

STK32A (NM_001287740) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	STK32A (NM_001287740) Human Untagged Clone
Tag:	Tag Free
Symbol:	STK32A
Synonyms:	YANK1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335918 representing NM_001287740. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGAGCGAACACTTCAAGAAAACCACAGTGTTTGATGAAAATGAAGATGTCAACTTTGACCCTTT
GAAATTTTGCAGCCATTGGGAAAGGCAGTTTTGGGAAGGTCTGCATTGTACAGAAGAATGATACCAAG
AAGATGTACGCAATGAAGTACATGAATAAACAAAAGTGCCTGGAGCGCAATGAAGTGAGAAATGTCTTC
AAGGAACCCAGATCATGCAGGGTCTGGAGCACCTTTCTGGTTAATTTGTGGTATTCTTCCAAGAT
GAGGAAGACATGTTTCATGGTGGTGGACCTCTGCTGGGTGGAGACCTGCCTTATCACCTGCAACAGAAC
GTCCACTTCAAGGAAGAAACAGTGAAGCTTTCATCTGTGAGCTGGTTCATGGCCCTGGACTACCTGCAG
AACCGAGCATCATTACAGGGATGAAGCCTGACAATATTTTACTTGACGAACATGGGCACGTGCAC
ATCACAGATTTCAACATTGCTGCGATGCTGCCAGGGAGACACAGATTACCACCATGGCTGGCACCAG
CCTTACATGGCACCTGAGATGTTTCAGCTCCAGAAAAGGAGCAGGCTATTCTTTGCTGTTGACTGGTGG
TCCCTGGGAGTGACGGCATATGAACTGCTGAGAGGCCGGAGACCGTATCATATTCGCTCCAGTACTTCC
AGCAAGGAAATTTGACACACGTTTGAGACGACTGTTGTAACCTTACCCTTTCGCTGGTTCACAGGAAATG
GTGTCACCTTCTAAAAAGCTACTCGAACCTAATCCAGACCAACGATTTTCTCAGTTATCTGATGTCCAG
AACTTCCCCTATATGAATGATATAAACTGGGATGCAGTTTTTTCAGAAAGAGGCTCATTCCAGGTTTCATT
CCTAATAAAGGCAGGCTGAATTTGATCCTACCTTTGAACTTGAGGAAATGATTTTGGAGTCCAAACCT
CTACATAAGAAAAAAGCGTCTGGCAAAGAAGGAGAAGGATATGAGGAAATGCGATTCTTCTCAGACA
TGCTTCTTCAAGAGCACCTTGACTCTGTCCAGAAGGAGTTCATAATTTTCAACAGAGAAAAAGTAAAC
AGGGACTTTAACAAAAGACAACCAATCTAGCCTTGAACAAAACCAAGACCCACAAGTGACAAATGGA
CAAATGGACACAGGACTCAGTGAGACTTTTTCAGACCTCGAAAGTTTCAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online >](#)

Plasmid Map:


ACCN: NM_001287740

Insert Size: 1224 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001287740.1](#)

RefSeq Size: 4923 bp

RefSeq ORF: 1224 bp

Locus ID: 202374

Cytogenetics: 5q32

Protein Families: Druggable Genome, Protein Kinase

MW: 47.5 kDa