

Product datasheet for **SC335932**

TEX13A (NM_001291277) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: TEX13A (NM_001291277) Human Untagged Clone
Tag: Tag Free
Symbol: TEX13A
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC335932 representing NM_001291277.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGGCCTTGAGACCTGAGGACCCAGTAGCGGGTTCGGCATAGCAACGTGGTGGCCTTCATCAACGAG
AAAATGGCCAGGCACACGAAAGGCCCGAGTCTATCTTGAGAATATATCCTTATCCTGGGAGAAGGTG
GAAGACAAGCTGAGGGCCATACTGGAGGACAGCGAGGTGCCAGTGAGGTCAAAGAGGCCTGCACCTGG
GGCAGCCTGGCCTTGGGAGTGCGCTTTGCCACAGGCAGGCACAGCTACAAAGGCACAGGGTGCGGTGG
CTGCACGGCTTCGCCAACTGCACAAATCAGCCGACAGGCCTTGGCATCAGACCTGAAGAAGCTCAGG
GAGCAGCAGGAGACGGAACGCAAGGAGGCGGCTCCCGGCTAAGAATGGCCAGACCGCTCGTGGAG
GTGCAGAAAGAGAGAGACAAGGAGCTGGTGTCTCCCATGAGTGGGAGCAGGGGGCAGGGTGCCAGGC
CTGGCCACTGCCGAGGGGTTTGACAGAAGGAGCAGCTGAGGAGGAAGAAGAGCGCGGTGGCTGCT
GCTGGTGTGCTGGAGGAAAAGGAGCAGAAGAAGAGCAGAGGGATGTGGAGGTTGTGGCTGCCCTGTG
GAGGCCATGGCTCCCCCTGTGGAGGCTGGGGCTGCCCCATGGAGACCCAGTTCCCCACGTGGAGGCC
AGGGCTGCCTCCATGGAGACCACAGAGAAGCTGGAGAGAATCCTCCTGCAGCTCCTTGGAGATGCTGAT
CAGGAAAAGTACACCTATTGGGGGCAGAAGGAGGAGATCTCCGGTCGGTCGAAACAGCCACATCTTAT
TTCTCTGGAACCACGAACCCCTGGTCCAGAGCCTCATCAGAACCTTCTCCTGTCCAGCTCCCTGCCTCA
TACTCATACTCATACTCAAGCCCTTTTCTCCTTCTCAGACATACCCACTATATCCCTCCACAAGCA
ACAGTCACAGCACCAGTTCCGCCTCAGCTGCCTTCCGACTGGGAGGCCTTGATACTAGCCTGTGGTCT
GATGGGGGGCCCCACAGAATAGACCATCAGGAGACCCAGAGACAGGAGATACTCCGAACCTCATCAG
CAAAGACCTCCAGTATATCGCAGGCCAGGGGACTGGGACTGCCCTTGGTGTAAACGCTGTGAATTTTCA
CGGAGGGATACTTGCTTCGACTGTGGGAAGGGAATCTGGCTGCAAAAACCTCAT

Restriction Sites: SgfI-MluI
ACCN: NM_001291277
Insert Size: 1230 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001291277.1</u>
RefSeq Size:	1383 bp
RefSeq ORF:	1230 bp
Locus ID:	56157
UniProt ID:	<u>Q9BXU3</u>
Cytogenetics:	Xq22.3
Protein Families:	Druggable Genome
MW:	45.6 kDa
Gene Summary:	This gene is similar to a mouse gene that is expressed in the testis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2014] Transcript Variant: This variant (2) differs in the 5' UTR, compared to variant 1. Both variants 1 and 2 encode the same protein.