

Product datasheet for **SC335935**

ZNF3 (NM_001278291) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ZNF3 (NM_001278291) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF3
Synonyms:	A8-51; HF.12; KOX25; PP838; Zfp113
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335935 representing NM_001278291. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTTGGCGGCTGCGCTCCTAAAGGCCAAGTCCCAGGAGCTGGTAACCTTTGAGGATGTAGCTGTGTAC
TTCATCCGGAAGGAGTGGAAAGCGTTTGGAACTGCTCAGAGGGACCTCTATAGAGATGTGATGCTGGAG
AATTACGGGAATGTGTTCTCACTGGATCGTGAGACCAGGACTGAAAATGATCAAGAAATTTCTGAAGAC
ACAAGATCACATGGGTCTACTGGGAAGATTTCAAAGGATATTTCTCAGGGTCTCAAGTTTAAAGAA
GCCTATGAACGAGAAGTCAGTCTGAAAAGGCCGCTGGGAACTCCCTGGAGAAAGACTGAACAGGAAA
ATGCCAGATTTTGGTCAAGTGACAGTTGAGGAGAAGCTAACCCAGGGGAGAGAGAAGCGAGAAATAT
AATGATTTTGGAAACAGCTTCACTGTGAATCCAACCTTATCTCACATCAGAGACTCCCGTGGGAGAC
AGACCCATAAGTGTGATGAATGTAGCAAGAGCTTAAATCGAACTTCAGACCTTATTCAACATCAGAGA
ATCCACACTGGGAAAAGCCCTATGAATGTAATGAGTGTGGAAAGCCCTTCAGCCAGAGCTCACACCTT
ATTCAGCATCAGAGAAATCCACACTGGGAAAAACCTTATGAATGTAGTGATTGTGGGAAAACCTTCAGC
TGTAGCTCTGCCCTATTCTGCATCGGAGGATCCACACGGGGAGAAACCTATGAATGTAATGAGTGT
GGGAAGACCTTCAGCTGGAGCTCCACCCTACCCACCATCAGAGAATCCACACTGGTGGAGAAACCTAC
GCCTGCAATGAATGTGGAAAGCCCTTCAGCAGGAGCTCAACCTTATTCCACATCAGAGAATCCACACT
GGAGAAAACCTATGAATGTAATGAATGTGGAAAGCCCTTCAGCCAGAGCTCACACCTCTACAGCAC
CAGAGAATCCACACTGGAGAGAAGCCCTACGAATGTATGGAATGTGGAGGAAAGTTTACCTACAGTTCA
GGCCTTATTACAGCATCAAAGAATCCACACCGGGGAGAACCCTATGAATGTAGTGTGGGAAAGCC
TTCAGGTACAGCTCGGCTCTTGTTCGCCATCAGAGAATTCACACTGGAGAGAAGCCTTTGAATGGGATC
GGCATGAGCAAAAGCTCCCTCAGAGTTACGACCGAGTTAAATATCAGAGAGTCCACGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online >](#)

Plasmid Map:


ACCN: NM_001278291

Insert Size: 1233 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278291.1](#)

RefSeq Size: 2605 bp

RefSeq ORF: 1233 bp

Locus ID: 7551

Cytogenetics: 7q22.1

Protein Families: Transcription Factors

MW: 47.1 kDa

Gene Summary:

Involved in cell differentiation and/or proliferation.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (6) differs in the 5' UTR, lacks the exon containing the translation start site, and initiates translation at an alternate start codon compared to variant 3. The resulting isoform (3) is shorter at the N-terminus compared to isoform 2. Variants 6, 7 and 13 encode the same isoform (3).