

Product datasheet for **SC335971**

MPP2 (NM_001278374) Human Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | MPP2 (NM_001278374) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | MPP2 |
| Synonyms: | DLG2 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >SC335971 representing NM_001278374. Blue=Insert sequence Red=Cloning site Green=Tag(s) |


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTGGGCATCCGAAGACAGCCGGAGAACATCTGGGTGTAACGTTCCGCGTGGAGGGCGGCGAGCTG
GTGATCGCGCGCATTCTGCATGGGGCATGGTGGCTCAACAAGGCTGCTGCATGTGGGTGACATCATC
AAGGAGGTGAACGGGCAGCCAGTGGGCAGTGACCCCGCGCACTGCAGGAGCTCCTGCGCAATGCCAGT
GGCAGTGTATCCTCAAGATCCTGCCAGCTACCAGGAGCCCATCTGCCCGCCAGGTATTTGTGAAA
TGTCACTTTGACTATGACCCGGCCGAGACAGCCTCATCCCTGCAAGGAAGCAGGCCTGCGCTTCAAC
GCCGGGACTTGCTCCAGATCGTAAACCAGGATGATGCCAACTGGTGGCAGGCATGCCATGTCGAAGGG
GGCAGTGTGGGCTATTCCAGCCAGCTGCTGGAGGAGAAGCGGAAAGCATTGTCAAGAGGGACCTG
GAGCTGACACAAACTCAGGGACCTATGCGGCAGCCTTTTCAGGAAAGAAAAGAAGCAATGATGTAT
TTGACCACCAAGAATGCAGAGTTTGACCGTCAAGCTGCTCATTATGAGGAGGTGGCCCGCATGCC
CCGTTCCGCCGAAAACCTGGTACTGATTGGGGCTCAGGGCGTGGGACGGCGCAGCCTGAAGAACAAG
CTCATCATGTGGGATCCAGATCGCTATGGCACCACGGTGCCTACACCTCCCGCGGGCCGAAAGACTCA
GAGCGGAAGGTCAGGGTTACAGCTTTGTGTCCCGTGGGAGATGGAGGCTGACGTCCGTGCTGGGGC
TACCTGGAGCATGGCGAATACGAGGGCAACCTGTATGGCACACGTATTGACTCCATCCGGGGCGTGGTC
GCTGCTGGGAAGGTGTGGTGTGGATGTCAACCCCGAGCGGTGAAGGTGCTACGAACGGCCGAGTTT
GTCCCTTACGTGGTGTTCATCGAGGCCAGACTTCGAGACCCTGCGGGCCATGAACAGGGCTGCGCTG
GAGAGTGAATATCCACCAAGCAGCTCACGGAGGCGGACCTGAGACGGACAGTGGAGGAGAGCAGCCGC
ATCCAGCGGGGCTACGGCACTACTTTGACCTCTGCCTGGTCAATAGCAACCTGGAGAGGACCTCCGC
GAGCTCCAGACAGCCATGGAGAAGCTACGGACAGACCCAGTGGGTGCCTGTCAGCTGGGTGACTGA
ACGGTACGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_001278374

Insert Size: 1242 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278374.1](#)

RefSeq Size: 4526 bp

RefSeq ORF: 1242 bp

Locus ID: 4355

UniProt ID: [Q14168](#)

Cytogenetics: 17q21.31

Protein Families: Druggable Genome

MW: 46.5 kDa

Gene Summary: Palmitoylated membrane protein 2 is a member of a family of membrane-associated proteins termed MAGUKs (membrane-associated guanylate kinase homologs). MAGUKs interact with the cytoskeleton and regulate cell proliferation, signaling pathways, and intracellular junctions. Palmitoylated membrane protein 2 contains a conserved sequence, called the SH3 (src homology 3) motif, found in several other proteins that associate with the cytoskeleton and are suspected to play important roles in signal transduction. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (6) lacks an in-frame coding exon and represents use of an alternate promoter, compared to variant 1. It encodes isoform 6 which has a shorter N-terminus compared to isoform 1.