

Product datasheet for **SC335974**

CCDC83 (NM_001286159) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC83 (NM_001286159) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC83
Synonyms:	CT148; HSD9; KP-CoT-23
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335974 representing NM_001286159. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAATACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAAACTCAGGAAAGCAAATAAAAAGGATACACATGACGGGCCACCAAAGAAATTAAGTGCCT
ACCAGTGAAGCACTTCTAGACTATCAATGTCAAATAAAGGAAGATGCCGTGGAGCAATTCATGTTTCAA
ATAAAGACACTTAGGAAAAAGAACCAAAAATATCATGAAAGAAATAGCCGCTTAAAGAAGAAGCAGATT
TGGCACATACGGCATCTACTAAAGGAAGTGAAGGAGAGGCGAGGGATTGCCAGTTGTAACAAGA
GAGGATGTTGAAGAAGCGATGAAGGAAAAATGGAAGTTTGAAGAGACCAGGAAAAAACTTGAGAGAT
ATGCGCATGCAAATAAGTAATGCTGAGAAACTATTTCTTGAGAACTCAGTGAAGGAATATTGGGAG
GAGTACAAGAATGTAGGGAGTGAACGACATGCTAACTCATTACCTCCTTACAAAATGACATCAACACA
GTTAAAGAGAATGCAGAGAAAATGTCAGAACACTATAAAATCACTCTGGAAGATACTAGAAAAGAAAATA
ATCAAGGAAACTTTGTTGCAACTGGACAAAAGGAATGGGCCACACAGAATGCTGTAAGCTCATT
GACAAGGGCAGTTATCTAGAGATCTGGGAGAAATGACTGGCTCAAAAAGAGGTTGCAATTCACAGGAAG
GAAGTTGAAGAATTAATAAATGCTATTCATGAACTGGAAGCAGAAAATTTGGTGCATTATTGATCAACTA
TCCAAGTGTAGACTTGTGGATCTCAAGATACCCAGGCGACTATATCTTACCCAAGCTGCTGGACTAGAA
GTGCCACCTGAAGAAATGCTTTTGAATTGCCAGAAACACATATAGAAGAGAAGTCAAGATTGCAACCC
ACAGAAGTAGAAAGTAGAGACTTGATGTCCTCATCAGATGAGAGCACTATCTTACATCTTAGTCATGAA
AATAGCATCGAAGATCTCCAGTATGTGAAGATAGATAAAGAGGAAAACAGGCACAGAGTTTGGGGAC
ACTGATATGAAGTACTTACTATATGAGGATGAGAAGGATTTCAAGGATTATGTAAGTGGGCCCCCTG
GGAGTGAAGCTTATGAGTGTGGAGAGCAAGAAAATGCCCATTCATTTTCAAGAGAAGGAAATCCAGTC
AAACTCTATAAAGATGTCAGGAGCCAGAAAGCCACATCAGATATAAGATGATGAAGTCTTTTCTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online >](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001286159

Insert Size: 1242 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286159.1](#)

RefSeq Size: 2276 bp

RefSeq ORF: 1242 bp

Locus ID: 220047

UniProt ID: [Q8IWF9](#)

Cytogenetics: 11q14.1-q14.2

MW: 48.9 kDa