

## Product datasheet for **SC335978**

### TUBA6 (TUBA1C) (NM\_001303117) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	TUBA6 (TUBA1C) (NM_001303117) Human Untagged Clone
Tag:	Tag Free
Symbol:	TUBA6
Synonyms:	bcm948; TUBA6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335978 representing NM_001303117. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAATACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCAAGTGACAAGACCATTGGGGGAGGAGATGATTCCCTCAACACCTTCTTCAGTGAACGGGTGCT
GGCAAGCATGTGCCCGGGCAGTGTGGTAGACTTGGAAACCCACAGTCATTGATGAAGTTCGCACTGGC
ACTTACC GCCAGCTCTTCCACCCTGAGCAACTCATCACAGGCAAGGAAGATGCTGCCAATAACTATGCC
CGAGGGCACTACACCATTGGCAAGGAGATCATTGACCTCGTGTGGACCGAATTCGCAAGCTGGCTGAC
CAGTGCACCGGTCTTCAGGGCTTCTGGTTTTCCACAGCTTTGGTGGGGAACTGGTTCTGGTTACC
TCGCTGCTCATGGAACGCTCTCAGTTGATTATGGCAAGAAGTCCAAGCTGGAGTTCTCCATTTACCCG
GGCCCCAGGTTTCCACAGCTGTAGTTGAGCCCTACAACCTCCATCCTCACCACCCACACCACCTGGAG
CACTCTGATTGTGCCTTCATGGTAGACAATGAGGCCATCTATGACATCTGTGCTAGAAAACCTCGATATC
GAGCGCCCAACCTACACTAACCTTAACCGCCTTATTAGCCAGATTGTGTCTCCATCACTGCTTCCCTG
AGATTTGATGGAGCCCTGAATGTTGACCTGACAGAATTCAGACCAACCTGGTGCCTACCCCGCATC
CACTTCCCTCTGGCCACATATGCCCTGTATCTCTGCTGAGAAAGCCTACCACGAACAGCTTACTGTA
GCAGAGATCACC AATGCTTGCTTTGAGCCAGCCAACCAGATGGTGAATGTGACCTCGCCATGGTAAA
TACATGGCTTGTGCCTGTTATACCGTGGTGACGTGGTTCCCAAAGATGCAATGCTGCCATTGCCACC
ATCAAAACCAAGCGTACCATCCAGTTTGTGGATTGGTGGCCCACTGGCTTCAAGTTGGCATTAAATTAC
CAGCTCCCACTGTGGTGCCTGGCGGAGACCTGGCCAAGGTACAGAGAGCTGTGTGCATGCTGAGCAAT
ACCACAGCTGTTGCCGAGGCCTGGGCTCGCCTGGACCACAAGTTTACCTGATGTATGCCAAGCGTGCC
TTTGTACTGTTACGTGGGTGAGGGGATGGAGGAAGGCGAGTTTTCAGAGGCCCGTGAGGACATGGCT
GCCCTTGAGAAGGATTATGAGGAGGTTGGAGCAGATAGTCTGACGGAGAGGATGAGGGTGAAGAGTAT
TAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001303117
<b>Insert Size:</b>	1245 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001303117.1</a></u>
<b>RefSeq Size:</b>	1825 bp
<b>RefSeq ORF:</b>	1245 bp
<b>Locus ID:</b>	84790
<b>Cytogenetics:</b>	12q13.12
<b>Protein Pathways:</b>	Gap junction, Pathogenic Escherichia coli infection
<b>MW:</b>	46.1 kDa
<b>Gene Summary:</b>	<p>Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (4) contains an alternate exon in the 5' end and initiates translation at an alternate start codon compared to variant 1. The encoded isoform (b) is shorter at the N-terminus compared to isoform a. Variants 2, 3, and 4 all encode the same isoform (b).</p>