

Product datasheet for **SC336024**

FLACC1 (NM_001289993) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: FLACC1 (NM_001289993) Human Untagged Clone
Tag: Tag Free
Symbol: FLACC1
Synonyms: ALS2CR12
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC336024 representing NM_001289993.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTACCCCAACCTCTCATCTACTGCACCTGCTGGGACCCCTGGAACCTGGGACCACGGAAGCTAATC
AAGACCCCTCACTACCACGCAAGAACTCCACAGGGAGTTCCAAGCTAACTCCTCTTGTACCAGCTCCA
AAAAATCACAATTACCTCCAACCAACAAAACCTGTTGTTTCCCCAAAAATGAAAATCCATTCAGCAAGG
CAAGAAGAGACTAATAAATCATTTTATGAAGTGATCAACGTGTCACCTGGCTATCAACTTGTTCCGAAT
CGGGAACAGATTTCTGTACCTTAGGGGATGAGATGTTTGATAGGAAAAAGCGGTGGGAATCGGAGATC
CCGGACAAAGGCAGATTTTCCAGGACCAACATCATTTCTGACCTAGAAGAGCAAATCTCAGAGCTGACA
GCAATAATTGAACAAATGAACAGAGACCACAGTCTGCCAGAAATTGCTCTCCAGTGAAATGGATCTC
CGCTGTGCTGAGATGAAACAGAACTTTGAAAACAAGAACAGGGAGCTCAAAGAGGCCCATGAAGCAGAA
CTCAGTGAGTTGGAGAACAACACTACAAAGCAGCCTTGAAGGCAGAGAAGTTGGCTGCCAAGAGAAGCTA
GAGGAGATGGGAAAAGAATACAAGTATTTGAAGAATATGTTTCGTACGTATCAGGACAGTATTTATGAT
GAAATGGAAGAGAAGTGGTCAAAACAGAAGGCGAAATGGAAGAAGGATGAGAAGTTCGAGCGAGAAAAAT
ATCCTGCTACAGCAAAAAAAAAAGATGACCAAAAAATTCGAAATGGAGTCAGGAGAAGAAGATAAGAAA
ATAAATGAATCCTGCAGTGCTGTCTTTGAGAACTTCATTCAAGAGAAGGAGGAGCTCTTGAACAACAT
CAAAGTGACACCTTGCAATTAGAAGAGCTGAGAAAAACCAAGAGGTGATGCAGGAAGAATTGCATGCA
CAAGCCCTTATCCTAGAGTCACTGAACACAAAACCTCTACTATACCCAGTTGGAACCTCCAGAAAGAGAAA
GCTATAGTGGGAAATCTGGAGAAAATGCTTCAAACCAAGTTTGCTGAAACTGAAGAAAAGTATAAGCAC
ACCATACAGATCCTGACGGAAGAGAACATTCATCTGAAGCAAAAGATAATTTCTAAGAAATGAAGAAATT
TGTGAAGGATGTTCTGGGAGATTGGCCTCTATTACTGTTTCTAAGGATGATTCTGACACTGTGCAAGAT
GGTAGCAAGAAAGGACAAGAATCATAA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001289993
Insert Size: 1269 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001289993.1</u>
RefSeq Size:	1785 bp
RefSeq ORF:	1269 bp
Locus ID:	130540
UniProt ID:	<u>Q96Q35</u>
Cytogenetics:	2q33.1
MW:	49.5 kDa