

## Product datasheet for **SC336041**

### ZNF385C (NM\_001242704) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF385C (NM_001242704) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF385C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336041 representing NM_001242704. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGTCTGGCTGCAGGCCCGGCCCTCCGGCGCCCCAGCCCCCTGCTGGCCTCCCTGCCCTGCCACC
CGGCCTCTGCAGCCCCGCTGGACTTCAAGCACTTGCTCGCCTTCCACTTCAATGGCGCTGCCCGCTC
AGTCTCTTCCCCAACTTCAGCACGATGGACCCGGTCCAGAAAGCTGTATCAGCCACACGTTTGGTGTC
CCCTCCCTCTGAAGAAGAAGCTGTTTCTTCTGTAACATCTGTACCTGAGGTTCAACTCAGCGAAC
CAGGCCGAGGCACATTATAAAGGCCACAACACGCCAGAAAACCAAGGCTGTCGAGGCTGCCAAGAGC
AAGCAGAGGCCACACACCCAGGCCAGGATGGGGCTGTAGTGTCCCAATCCCAACGCTGGCCAGTGGA
GCCCTTGAGAGCCACAGAGTAAAGTTCCTGCAGCCCTCCTCTTGGGCCCTCCACTCCAGCCACCACCA
ACTCCAGACCCTACATGCAGGGAGCCAGCCACTCAGAGCTCTTGATGCTGCCTCCTCATCCTTCTCT
TCTCCTGCCACCTTGCTCCCCAGAGCCTGGGAGAGAGGCACCGGGCCTGAGCCAGCGGCAGCTGCC
GTGGGAAGCAGCATGAGTGGGAAGGCAGGAGTGAGAAGGGGCACCTCTACTGCCACGTGTAAGGTG
ACAGTGAACCTCGGCCTCCAGCTTCAGGCTCACAACACAGGAGCCAAGCACCGGTGGATGATGGAAGGT
CAGCGAGGGGCTCCCCGAGGAGCCGGGGCCGCCGGTGTCCAGGGGAGGTGCCGACACAAAGCCAAG
AGAGTACAGGGGGCCGGGGCGGCCAGGGGCCAGCCCTGCCTTCCACTGTGCTCTCTGTAGCTC
CAGGTCAATTCAGAGACCAACTGAAGCAGCACATGAGCAGCAGGAGGCACAAAGACCGCCTGGCCGGG
AAGACCCCAAGCCCTCCAGCCAGCACAGCAAGCTGCAGAAGCACGCAGCGCTGGCTGTGAGTATCCTC
AAGTCTAAACTGGCCTTGAGAAGCAACTACCAAGACGTTGGCAGCCCGCTTCTGCCAGCCCGCTC
CCCACCGCAGCCACTGCCATCTGTGCTCTGCCAGGGCCCTGGCCCTCCGCCCTGCCCTACAGCAGCC
ACTACCCTCTTCCCGCTCCCATCTGGGCCAGCTCTGTTTCGCACCCAGCAGGAGCTGTCCGCCCT
GCCACAGGACCTATCGTCTTGGCCCTTAG
ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online >](#)

<b>ACCN:</b>	NM_001242704
<b>Insert Size:</b>	1275 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001242704.1</a></u>
<b>RefSeq Size:</b>	2596 bp
<b>RefSeq ORF:</b>	1275 bp
<b>Locus ID:</b>	201181
<b>Cytogenetics:</b>	17q21.2
<b>MW:</b>	44.2 kDa