

Product datasheet for SC336051

MEK7 (MAP2K7) (NM_001297556) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MEK7 (MAP2K7) (NM_001297556) Human Untagged Clone
Tag:	Tag Free
Symbol:	MAP2K7
Synonyms:	JNKK2; MAPKK7; MEK; MEK 7; MKK7; PRKMK7; SAPKK-4; SAPKK4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336051 representing NM_001297556. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGCGCTCCTCCCTGGAACAGAAGCTGTCCCCTGGAAGCAAAGCTGAAGCAGGAGAACCGGGAG
GCCCGCGGAGGATCGACCTCAACCTGGATATCAGCCCCAGCGGCCAGGCCACCCCTGCAGCTCCCG
CTGGCCAACGATGGGGCAGCCGCTCGCCATCCTCAGAGAGCTCCCCGAGCAGCCACGCCCCCGCC
CGCCCCGCCACATGCTGGGGCTCCCGTCAACCCTGTTACACCCCGCAGCATGGAGAGCATTGAGATT
GACCAGAAGCTGCAGGAGATCATGAAGCAGACGGGCTACCTGACCATCGGGGCCAGCGCTACCAGGCA
GAAATCAACGACCTGGAGAACTGGGCGAGATGGGCAGCGCACCTGCGGCCAGGTGTGGAAGATGCGC
TTCCGGAAGACCGGCCAGTTCATTGCCGTTAAGCAAATGCGGCGCTCCGGGAACAAGGAGGAGAACAAG
CGCATCCTCATGGACCTGGATGTGGTGCTGAAGAGCCACGACTGCCCTACATCGTGCAGTGCTTTGGG
ACGTTTCATCACCAACACGGACGCTTTCATCGCCATGGAGCTCATGGGCACCTGCGCTGAGAAGCTCAAG
AAGCGGATGCAGGGCCCCATCCCCGAGCGCATTCTGGGCAAGATGACAGTGGCGATTGTGAAGGCGCTG
TACTACCTGAAGGAGAAGCAGGTGTATCCACCGCAGCTCAAGCCCTCCAACATCCTGCTGGACGAG
CGGGGCCAGATCAAGCTTGCAGCTTCCGATCAGCGCCCGCTGGTGGACTCCAAGCCAAGACCGCGG
AGCGCCGGCTGTGCCCTACATGGCACCCGAGCGCATTGACCCCCAGACCCCAAGCCGGACTAT
GACATCCGGGCGACGTATGGAGCCTGGGCATCTCGTTGCCCTGCCCTCCCAGGTGGAGCTG
GCAACAGGACAGTTCCCTACAAGAACTGCAAGACGAGCTTTGAGGTCTCACAAAGTCTACAGGAA
GAGCCCCCGCTTCTGCCCGACACATGGCTTCTCGGGGACTTCCAGTCTTCTGCAAAGACTGCCTT
ACTAAAGATCACAGGAAGAGACCAAGTATAATAAGCTACTTGAACACAGCTTCAAGCGCTACGAG
ACGCTGGAGGTGGACGTGGCGTCTGTTCAAGGATGTCATGGCGAAGACTGAGTACCAGCGGACTAGC
GGCGTCTGAGCCAGCCCCACCTGCCCTTCTCAGGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001297556

Insert Size: 1281 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001297556.1](#)

RefSeq Size: 3498 bp

RefSeq ORF: 1281 bp

Locus ID: 5609

UniProt ID: [O14733](#)

Cytogenetics: 19p13.2

Protein Families: Druggable Genome, Protein Kinase

Protein Pathways:	ErbB signaling pathway, Fc epsilon RI signaling pathway, GnRH signaling pathway, MAPK signaling pathway, Neurotrophin signaling pathway, T cell receptor signaling pathway, Toll-like receptor signaling pathway
MW:	48.2 kDa
Gene Summary:	<p>The protein encoded by this gene is a dual specificity protein kinase that belongs to the MAP kinase kinase family. This kinase specifically activates MAPK8/JNK1 and MAPK9/JNK2, and this kinase itself is phosphorylated and activated by MAP kinase kinase kinases including MAP3K1/MEKK1, MAP3K2/MEKK2, MAP3K3/MEKK5, and MAP4K2/GCK. This kinase is involved in the signal transduction mediating the cell responses to proinflammatory cytokines, and environmental stresses. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]</p> <p>Transcript Variant: This variant (2) lacks an alternate in-frame exon in the 5' coding region and uses an alternate in-frame splice site in the 3' coding region, compared to variant 1. It encodes isoform 2, which is shorter than isoform 1.</p>