

## Product datasheet for **SC336052**

### TTC39A (NM\_001297662) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	TTC39A (NM_001297662) Human Untagged Clone
Tag:	Tag Free
Symbol:	TTC39A
Synonyms:	C1orf34; DEME-6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336052 representing NM_001297662. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACCTCGGCTGGCGCGCCCCAGGAGCCCTGCCCGCGGGGACTCCTGAGAGCAGCCTCCATGAGGCC
CTGGACCAGTGCATGACCGCCCTGGACCTTCTCCTACCAACCAGTTCTCAGAAGCACTCAGCTACCTC
AAGCCCAGAACCAAGGAAAGCATGTACCACTCACTGACATATGCCACCATCCTGGAGATGCAGGCCATG
ATGACCTTGGACCTCAGGACATCTGCTTCCCGCAACATGATGAAGGAGGCACAGATGCTGTGTCAG
AGGCACCGGAGGAAGTCTTCTGTAAACAGATTCTTCAGCAGCCTGGTGAACCGCCACGCTGGGCCAA
TTCAGTGAAGAGGAAATCCACGCTGAGGTCTGCTATGCAGAGTGCCTGCTGCAGCGAGCAGCCCTGACC
TTCTGCAGGGTTCTCACACGGAGGGGAGTCAGGCCAGAGCCTTGATGATCCCTCTCACGCCCTGC
AGCTGCCACCTGGGCCAGGCCGTCAGCATCTTTTCTCCTGCAGGACGAGAATCGGTGAGCTTCAAC
AAAGGCGGCATCAAAGTTCGAAACAGCTACCAGACCTACAAGGAGCTGGACAGCCTTGTTTCACTCTCA
CAATACTGCAAGGGTGAGAACCCCGCACTTTGAAGGAGGAGTGAAGCTTGGTGTAGGGGCCCTTCAAC
CTGACACTGTCCATGCTTCTACTAGGATCCTGAGGCTGTTGGAGTTTGTGGGTTTTTTCAGGAAACAAG
GACTATGGGCTGCTGCAGCTGGAGGAGGGAGCGTCAGGGCACAGCTTCCGCTCTGTGCTCTGTGTCATG
CTCCTGCTGTGCTACCACACCTTCTCACCTTCGTGCTCGGTACTGGGAACGTCAACATCGAGGAGGCC
GAGAAGCTTTGAAGCCCTACCTGAACCGGTACCCTAAGGGTGCCATCTTCTGTTTTCAGGGGAGG
ATTGAAGTCATTAAGGCAACATTGATGCAGCCATCCGGCGTTTTTCAGGAGTGTGTGAGGCCAGCAG
CACTGGAAGCAGTTCACACATGTGCTACTGGGAGCTGATGTGGTCTTACCTACAAGGGCCAGTGG
AAGATGCTCACTTCTACGCCACTGCTCAGCAAGGAGAACTGCTGGTCCAAGGCCACCTACATTTAC
ATGAAGGCCGCTACCTCAGCATGTTTGGGAAGGAGGACCACAAGCCGTTCCGGGACGACGAAGTGGAA
TTATTTCCGGTGGTCCCACATCGTCTCCCACTGTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: Sgfl-Mlul

**Plasmid Map:**



ACCN: NM\_001297662

Insert Size: 1281 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001297662.1](#)

RefSeq Size: 1892 bp

RefSeq ORF: 1281 bp

Locus ID: 22996

Cytogenetics: 1p32.3

MW: 48.3 kDa