

Product datasheet for SC336056

RUSC1 (NM_001278228) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RUSC1 (NM_001278228) Human Untagged Clone
Tag:	Tag Free
Symbol:	RUSC1
Synonyms:	NESCA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336056 representing NM_001278228. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGGAGAGTAATGGAGCTCCGCGGGGGTGCAGGAATGGTTGGCGGGAGGTGGCTGGGAGGTGTCTG
GAGGGATCTGGAGCCATCGGCATCGCGCCACCTCCAGGCCAGTTGGGTGATAGCCGGCTGAGCCCG
GATGTGGGGCACCTGGTGTGACCACCTCTGCCCGCCCTCCACGCCCTGGTGGCGGACGGGCTGAAG
CCTTCCGGAAGGACCTCATCCCGGGCAGCGCAGGAGCAGCCCTGGAGCGTGGTGGAGGCGTCGGTG
AAGCCAGGCTCCAGCACCCTCCCTTGAACCTGTATAGCCAGGTCAGCCGTAGCCCCGCTGAGC
AGCAGCCGTAGCCGCTTCCATGCCTTATCCTGGGCTCCTCAACACCAAGCAGTTGGAGCTGTGGTTT
TCCAGTCTCCAGGAAGATGCAGGCTGCTCCTCCTGTACCTGCCAACAGGATTTTCTCCCTGGCC
CGCGGTGGTTGTCCCTCCCTGTCCACAGAGCTGCTGCTCCTGCTGCAGCCATTGTCGGTGTCACTTTC
CACCTGGACCTGCTCTTTGAGCACCACCACCACCTGCCCTGGGCCACCTCAGGCCCTGCCCTCCA
GGCCACCTCCAGCTCTGCAGCAGACTATGCAAGCCATGCTGCACCTTTGGGGCCGGCTGGCCAGAGC
CTTCGGGGGACTTCCAAGGAAGCTGCTTCCAGACCCTCTGACTCTCAAACCTTCCCACACCAGGGAGC
TGGTGGGAGCAGTTGACCAGGCTCCCGGGTCTATGCCTCTGGGGGCACTGAGGGCTTTCCTCTTCC
CGATGGGCACCGGGCGTCATGGGACTGCAGCTGAAGAAGGTGCACAGGAGAGACCCCTGCCACAGAT
GAGATGGCACCAGGCAGGGCCCTCTGGTTGGGAAGACTATTTGGAGTGCCTGGGGCCCGCAGAAAAAT
GAGAATGGAGCCCTAAAGTCCAGGAGACCATCTAGCTGGCTGCCCGGACAGTGTGTGGCTCTT
GTGAAGCGGGGGGCACCTCCCGAGATGCCTTCTCCTCAGGAGCTTGGGGCTCAGCACCCAGGATGGTG
CAAACCCATAGGGCAGTGCAGGCTCTGTGTGATCACACTGCTGCAAGACCTGACCAGTTGAGCTCCGG
CGTGGGGAAGTGTGCTGCTCATCACCACAGTGGATGAGGACTGGCTCCGCTGTGGGCGGGATGGCATG
GAGGCTGTGGTGTGGGTATACCTCCCTTGTCTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online >](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001278228

Insert Size: 1284 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278228.1](#)

RefSeq Size: 2705 bp

RefSeq ORF: 1284 bp

Locus ID: 23623

UniProt ID: [Q9BVN2](#)

Cytogenetics: 1q22

MW: 46.1 kDa

Gene Summary:

Putative signaling adapter which may play a role in neuronal differentiation. May be involved in regulation of NGF-dependent neurite outgrowth. Proposed to play a role in neuronal vesicular trafficking, specifically involving pre-synaptic membrane proteins. Seems to be involved in signaling pathways that are regulated by the prolonged activation of MAPK. Can regulate the polyubiquitination of IKBKG and thus may be involved in regulation of the NF-kappa-B pathway.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (6) differs in the 5' UTR and contains multiple differences in the coding region, compared to variant 1. It initiates translation at an alternate start codon. The encoded isoform (f) is shorter and has a distinct N-terminus, compared to isoform a.