

Product datasheet for **SC336062**

ELP3 (NM_001284224) Human Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | ELP3 (NM_001284224) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | ELP3 |
| Synonyms: | KAT9 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >SC336062 representing NM_001284224. Blue=Insert sequence Red=Cloning site Green=Tag(s) |


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCCACACATCAGTTTTACAGGAAATATATGTGTCCAACCTCCATGAGAGCTATCCGTGCCAGATAT
GACCCCTTCTACAGACAAGACACCGAATAGAACAGTTAAACAACCTGGTCATAGTGTGGATAAAGTG
GAGTTTATTGTGATGGGTGGAACGTTTATGGCCCTTCCAGAAGAATACAGAGATTATTTTATTGAAAT
TTACATGATGCCTTATCAGGACATACTCCAACAATTTACGAGGCAGTCAAGTATTCTGAGAGAAGC
CTCACAAAGTGTATTGGAATTACTATTGAAACCAGACCAGATTACTGCATGAAGCGACATTTAAGTGAC
ATGTTGACCTATGGCTGCACAAGGCTGGAGATTGGGGTGCAGAGTGTATGAAGATGTGGCTAGAGAC
ACCAACAGGGGCCACACTGTGAAGGCAGTGTGTGAGTCATTTACCTGGCCAAAGATTCCGGTTTTAAA
GTGGTGGCCCATATGATGCCTGACCTGCCAAACGTTGGGACTAGAAAGAGACATTGAACAGTTCACAGAG
TTTTTTGAGAACCCTGCTTTTCGTCCCGATGGGCTGAAACTCTATCCTACCCTGGTGATTCGTGGGACC
GGGCTTTATGAGCTTTGGAAATCAGGAAGATAAAGAGTTACTCTCCTAGTGACCTGGTTGAATTGGTG
GCTCGGATCCTAGCCCTCGTGCCCTCATGGACTCGAGTGTACCGAGTACAGAGGGATTTCCAATGCCT
TTAGTTAGCTCAGGAGTAGAGCATGGTAACCTGAGAGAGCTGGCACTTGAAGAATGAAAGACCTCGGA
ATACAGTGTGAGATGTGAGAACCAGAGAAGTTGGAATCCAAGAAATTCATCACAAAGTACGGCCATAC
CAGGTTGAATTGGTAAGGAGAGATTATGTTGCAAATGGTGGCTGGGAAACATTTGTGCATACGAAGAC
CCAGATCAAGACATTTTATTGGCTCCTACGATTACGCAAGTGTTCAGAAGAAACTTTCCGTTTTCGAA
TTGGGTGGAGGTGTCTCCATAGTACGAGAGCTGCATGTGTATGGGAGTGTGGTCCCTGTGAGCAGCCGG
GATCCTACTAAATTTAGCATCAGGGATTTGGCATGCTGCTGATGGAGGAAGCAGAAAGAATAGCTAGA
GAAGAACATGGGTCTGGGAAAAATCGCTGTGATATCAGGGGTGGCACCAGGAATTATTATAGAAAGATC
GGCTACAGATTACAAGGCCCGTACATGGTGAAGATGCTGAAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001284224

Insert Size: 1287 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001284224.1](#)

RefSeq Size: 3032 bp

RefSeq ORF: 1287 bp

Locus ID: 55140

UniProt ID: [Q9H9T3](#)

Cytogenetics: 8p21.1

MW: 49.4 kDa

Gene Summary:

ELP3 is the catalytic subunit of the histone acetyltransferase elongator complex, which contributes to transcript elongation and also regulates the maturation of projection neurons (Creppe et al., 2009 [PubMed 19185337]).[supplied by OMIM, Apr 2009]

Transcript Variant: This variant (4) differs in the 5' UTR and lacks a portion of the 5' coding region, compared to variant 1. These differences cause translation initiation at a downstream AUG and result in an isoform (4) with a shorter N-terminus, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments. Variants 4 and 5 encode the same isoform (5).