

Product datasheet for **SC336067**

MRPS27 (NM_001286748) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MRPS27 (NM_001286748) Human Untagged Clone
Tag:	Tag Free
Symbol:	MRPS27
Synonyms:	MRP-S27; S27mt
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336067 representing NM_001286748. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTGCCTCCATAGTGCGGCGCGGGATGCTCCTGGCGCGCAAGTGGTTCTTCCTCAGCTCTCTCCT
GCAGGTAAGAGATACCTGCTTTCTTCAGCCTATGTAGACAGCCACAAATGGGAAGCAAGAGAAAAAGAA
CATTACTGTCTTGCTGATCTTGCATCTTTAATGGATAAAACATTTGAGAGAAAGTTGCCTGTTAGTTCT
TTAACAAATACACGGCTTATAGACAACATTTCTCTCGGGAAGAGATAGATCATGCAGAGTATTACCTT
TACAAAAAGGGAGAATGCTCATCCAGCAATCCACAGAATTATTCCAGTTTCGACACAGCCCAACTGC
TGGTACCTGAGAACTGGACTATCCACACCTGGATTAGGCAGTGTCTAAAATATGATGCACAAGACAAA
GCCCTATATACCTTGTAAATAAGTTCAATATGGAAATTTTCCAGATAACTTTACATTTCAATTTACTG
ATGGATTCTTTCATAAAGAAAGAAAATTACAAAGATGCTTTATCTGTGGTTTTTGAGGTCATGATGCAA
GAAGCCTTTGAAGTGCCTTCCACCCAATTCTCTCCCTCTATGTTTTATTTTCATTGCCTGGCAAAGAAG
ACAGACTTCAGTTGGGAAGAGGAGAGAACTTTGGTGCATCCCTTTTGTCTCCAGGCCTAAAAACAAAAG
AACTCAGTGGGTTTCAGTCCCAGTTGTATGGCTATGCACTTCTTGGGAAGGTGGAGTTGCAGCAAGGG
CTACGGGCTGTGTACCACAACATGCCTCTGATATGGAAACCAGGCTACCTTGACAGAGCCCTTCAAGTG
ATGGAGAAAGTGGCTGCCTCCCCAGAAGACATAAAGCTGTGTAGAGAAGCGCTCGATGTGCTGGGTGCA
GTGCTGAAGGCTCTGACTTCAGTTCAGTGTAGGGCTTCAGAGGAGCAGTCCCAAAATGATGAAGACAACCAG
GGGTGAGAAAAACTGGTGGAGCAGTTAGACATCGAGGAAACAGAGCAGTCCAAGCTTCTCAATACCTG
GAACGATTTAAGGCCTTACATTCTAAGCTTCAAGCTCTGGGCAAAAATTGAGTCAGAAGGTCTTTTAAAGT
CTGACCACCCAGCTTGTCAAGGAAAACTCTCCACCTGTGAAGCAGAGGACATCGCCACCTATGAGCAG
AATCTGCAGCAGTGGCATCTAGACCTTGTACAGTTGATCCAGAGAGAACAGCAACAGAGGGAGCAAGCG
AAGCAGGAGTACCAGGCTCAGAAAGCAGCAAAGGCATCTGCCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



View online »

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_001286748

Insert Size: 1287 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286748.1](#)

RefSeq Size: 2839 bp

RefSeq ORF: 1287 bp

Locus ID: 23107

UniProt ID: [Q92552](#)

Cytogenetics: 5q13.2

MW: 49.1 kDa

Gene Summary:

Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein that may be a functional partner of the death associated protein 3 (DAP3). Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Nov 2013]

Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.