

Product datasheet for **SC336097**

CD36 (NM_001289908) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CD36 (NM_001289908) Human Untagged Clone
Tag:	Tag Free
Symbol:	CD36
Synonyms:	BDPLT10; CHDS7; FAT; GP3B; GP4; GPIV; PASIV; SCARB3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336097 representing NM_001289908. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGCTGTGACCGAACTGTGGCTCATCGCTGGGGCTGTCATTGGTGCTGTCCTGGCTGTGTTTGGG
GGTATTCTAATGCCAGTTGGAGACCTGCTTATCCAGAAGACAATAAAAAGCAAGTTGTCCTCGAAGAA
GGTACAATTGCTTTTAAAAATTGGGTTAAACAGGCACAGAAAGTTACAGACAGTTTTGGATCTTTGAT
GTGCAAAATCCACAGGAAGTGATGATGAACAGCAGCAACATTCAAGTTAAGCAAAGAGGTCCTTATACG
TACAGAGTTTCGTTTTCTAGCCAAGGAAAATGTAACCCAGGACGCTGAGGACAACACAGTCTCTTTCTG
CAGCCCAATGGTGCCATCTTCGAACCTTACTATCAGTTGGAACAGAGGCTGACAACCTTACAGTTCTC
AATCTGGCTGTGGCAGCTGCATCCCATATCTATCAAAATCAATTTGTTCAAATGATCCTCAATTCATT
ATTAACAAGTCAAAATCTTCTATGTTCCAAGTCAGAACTTTGAGAGAAGTGTATGGGGCTATAGGGAT
CCATTTTTGAGTTTGGTCCGTACCCTGTTACTACCACAGTTGGTCTGTTTTATCCTTACAACAATACT
GCAGATGGAGTTTATAAAGTTTTCAATGGAAAAGATAACATAAGTAAAGTTGCCATAATCGACACATAT
AAAGGTAAGGTCATCTATGCTGATTTGAATCCGACGTTAATCTGAAAGGAATCCCTGTGTATAGA
TTTGTCTTCCATCCAAGGCCTTGCCTCTCCAGTTGAAAACCCAGACAACACTATTGTTTCTGCACAGAA
AAAATATCTCAAAAAATTGTACATCATATGGTGTGCTAGACATCAGCAAATGCAAAGAAGGGAGACCT
GTGTACATTTCACTTCTCATTCTTGTGATGCAAGTCTGATGTTTCAGAACCTATTGATGGATTAAC
CCAAATGAAGAAGAACATAGGACATACTTGGATATTGAACCTATAAAGTGGATTCACTTTACAATTTGCA
AAACGGCTGCAGGTCAACCTATTGGTCAAGCCATCAGAAAAAATCAAGTATTAAGAATCTGAAGAGG
AACTATATTGTGCCTATTCTTTGGCTTAATGAGACTGGGACCATTTGGTGTGAGAAGGCAAAACATGTTT
AGAAGTCAAGTAACTGGAAAAATAAACCTCCTTGGCCTGATAGAAATGATCTTACTCAGTGTGGTGTG
GTGATGTTTGTGCTTTTATGATTTTCAATTGTGCATGCAGATCGAAAACAATAAAATTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: Sgfl-Mlul

Plasmid Map:



ACCN: NM_001289908

Insert Size: 1302 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001289908.1](#)

RefSeq Size: 1665 bp

RefSeq ORF: 1302 bp

Locus ID: 948

UniProt ID: [P16671](#)

Cytogenetics: 7q21.11

Protein Families: Druggable Genome, Transmembrane

Protein Pathways:	Adipocytokine signaling pathway, ECM-receptor interaction, Hematopoietic cell lineage, PPAR signaling pathway
MW:	48.6 kDa
Gene Summary:	<p>The protein encoded by this gene is the fourth major glycoprotein of the platelet surface and serves as a receptor for thrombospondin in platelets and various cell lines. Since thrombospondins are widely distributed proteins involved in a variety of adhesive processes, this protein may have important functions as a cell adhesion molecule. It binds to collagen, thrombospondin, anionic phospholipids and oxidized LDL. It directly mediates cytoadherence of <i>Plasmodium falciparum</i> parasitized erythrocytes and it binds long chain fatty acids and may function in the transport and/or as a regulator of fatty acid transport. Mutations in this gene cause platelet glycoprotein deficiency. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Feb 2014]</p> <p>Transcript Variant: This variant (6) lacks two consecutive exons in the 3' coding region, compared to variant 3. The resulting isoform (2) lacks an internal segment, compared to isoform 1.</p>