

Product datasheet for SC336108

FAM222B (NM_001288639) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM222B (NM_001288639) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM222B
Synonyms:	C17orf63
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336108 representing NM_001288639. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGAACCCCCAGTGGCACCCCTATGCTACTGTGGCACCCAGCACTTTAGCCCACCCAGGCCAGGCT
CTGGCCCCCAGCAGGCCCTGCAGCATGCACAGACCCCTGGCCCATGCCCTCCCCAGACGCTGCAGCAC
CCTCAGGGTATCCCGCCACCCAGGCAGTGTCCCACCCTCAGAGCCTCCAGCAGCCTCAGGGCCTGGGC
CACCCCTAGCCCATGGCCAAACCCAGGGCTTGGTCCACCCTCAGGCCCTGGCTCACCAGGGTCTCCAG
CACCCCAACAATCCCTTGTGCATGGAGCCGGAAGATGCCAGACTCAGATGCCCCCGAATGTGACC
GTGTCTACCTCAACTATCCCTTTCAATGGCGGCCACTCTGCAGCACAGCCAGCCTCCGGACCTGAGT
AGCATCGTGCACCAGATCAACCAGTTTTGCCAGACGAGGGCAGGCATCAGCACTACCTCAGTGTGTGAG
GGCCAGATCGCCAACCCAGCCCATAGTCGAGTCTGCTCATCAATGCAAGCACCCGGGTGTCGACC
CACAGCGTCCCCACACCAATGCCTTTCATGTGTGGTCAATCCCATGGAGCACACCCACGCGGCCACCGCC
CGGTTGCTGTGTCAGGTCTGTCAACCTGCCACAGGCATCTCTCGCGTCCCCACTGGCTACCCTAGC
GACCTCAAGCCAGTCACTGGAACCAGCACCAGCTGGCCACCTACAGCAGATGTGCAGCGAGGCTAGT
GGGACGCCAGCCCTGGCCTGACAGGCAAGCATGCGGCAGGACCGGAGTTGGCAGGGCCTGGCTTTGTG
GGCAAGGCCCTGCCTACCCGAGGAACCTGCCTGGCGCAGTCTTCCATCTGAAGCCACCCCTGGAA
AAGCCGACCCCATCCCCACCAAGTCAACGGCATGGCAGCCCATTTGGCCTACCCCAATGGTCACTACTTC
CAACCCCTGTGGAACAACATTCTGCCCACTCCCAATAGCGACAGCTCGGGGTCTCAGGACCTTGCCATG
CCGTTCCACGGTGGGCAGCCACAGGTGACCCCTCGACTGTGCGGCAGTCCCGGGGCCACTACCGA
GCAGGGACCGGGGGCGTCCAGTGGCAAGCCAGAACAGCTTGATGCAAACAGTGGATTACCTAAGTGGG
GATTTCCAACAGGCCTGCTTCCGAGAACAGAGCCTGGCCATGCTGAGCAAGGCCACCGAGCCCTGGC
AACCGAGCCCCGATCCACAGAGAGTCGAAGTCTTCATATTCAGCACCCAGGGTATAGATAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: Sgfl-Mlul

Plasmid Map:



ACCN: NM_001288639

Insert Size: 1305 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001288639.1](#)

RefSeq Size: 4101 bp

RefSeq ORF: 1305 bp

Locus ID: 55731

UniProt ID: [Q8WU58](#)

Cytogenetics: 17q11.2

MW: 45.6 kDa