

## Product datasheet for **SC336115**

### AAMP (NM\_001302545) Human Untagged Clone


#### Product data:

Product Type:	Expression Plasmids
Product Name:	AAMP (NM_001302545) Human Untagged Clone
Tag:	Tag Free
Symbol:	AAMP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336115 representing NM_001302545. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGTCCGAATCGAAAGCGGGGCTGCTGCTGACACCCCCCACTGGAGACCCTAAGCTTCCATGGT
GATGAAGAGATTATCGAGGTGGTAGAATTGATCCCGTCCGCCGACCCAGCAGATGACCTGGCCAG
GAGATGGAAGATGTGGACTTTGAGGAAGAAGAGGAGGAAGAGGGCAACGAAGAGGGCTGGTTCTAGAA
CCCCAGGAAGGGGTGGTCGGCAGCATGGAGGGCCCCGACGATAGCGAGGTCACCTTTCATTGCACTCA
GCATCTGTGTTTTGTGTGAGCCTGGACCCCAAGACCAATACCTTGGCAGTGACCGGGGTGAAGATGAC
AAAGCCTTCGTATGGCGGCTCAGCGATGGGGAGCTGCTTTGAGTGTGCAGGCCATAAAGACTGTGTG
ACTTGTGCTGGTTTCAGCCATGACTCCACTCTAGTGGCCACAGGGGACATGAGTGGCCTCTTGAAAGTG
TGGCAGGTGGACACTAAGGAGGAGTCTGGTCCTTTGAAGCGGGAGACCTGGAGTGGATGGAGTGGCAT
CCTCGGGCACCTGTCTGTTGGCGGGCACAGCTGACGGCAACACCTGGATGTGGAAGTCCCGAATGGT
GACTGCAAGACCTTCCAGGGTCCCAACTGCCAGCCACCTGTGGCCGAGTCTCCCTGATGGGAAGAGA
GCTGTGGTAGGCTATGAAGATGGGACCATCAGGATTTGGGACCTGAAGCAGGGAAGCCCTATCCATGTA
CTGAAAGGGACTGAGGGTACCAGGGCCCACTCACCTGTGTTGCTGCCAACCAGGATGGCAGCTTGATC
CTAACTGGCTCTGTGGACTGCCAGGCCAAGCTGGTCAGTGCCACCACCGCAAGGTGGTGGTGTTTTT
AGACCTGAGACTGTGGCCTCCCAGCCCAGCCTGGGAGAAGGGGAGGAGAGTGAGTCCAACCTGGTGGAG
TCCTTGGGCTTCTGCAGTGTGATGCCCTGGCAGCTGTTGGCTACCTGGATGGGACCTTGGCCATCTAT
GACCTGGCTACGCAGACTCTTAGGCATCAGTGTGACACCAGTCCGGCATCGTGCAGCTGCTGTGGGAG
GCAGGCACTGCCGTGGTATATACCTGCAGCCTGGATGGCATCGTGCCTCTGGGACGCCGGACCGGC
CGCCTGCTTACTGACTACCGGGGCCACAGCGCTGAGATCCTGGACTTTGCCCTCAGCAAAGATGCCTCC
CTGGTGGTGACCACGTGAGGAGACCACAAAGCGAAAGTATTTTGTGTCCAAGGCCTGACCGTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



**Plasmid Map:**


**ACCN:** NM\_001302545

**Insert Size:** 1308 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001302545.1](#)

**RefSeq Size:** 1845 bp

**RefSeq ORF:** 1308 bp

**Locus ID:** 14

**Cytogenetics:** 2q35

**MW:** 46.8 kDa

**Gene Summary:**

The gene is a member of the immunoglobulin superfamily. The encoded protein is associated with angiogenesis, with potential roles in endothelial tube formation and the migration of endothelial cells. It may also regulate smooth muscle cell migration via the RhoA pathway. The encoded protein can bind to heparin and may mediate heparin-sensitive cell adhesion.

[provided by RefSeq, Oct 2014]

Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).