

Product datasheet for SC336116

ZNF285 (NM_001291490) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF285 (NM_001291490) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF285
Synonyms:	ZNF285A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336116 representing NM_001291490. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTGAAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACCGAGCCCCAGAACTCTCAGGGAAGATATAAGGGAATTTACATGGAAGAGAAATTGTACAGACGT
GCTCAGCATGATGACAGCCTCAGTTGGACCTCATGTGATCATCATGAGTCCCAAGAAATGTAAGGAGAG
GACCCTGGTAGACATCCCAGCTGTGGGAAAACTTGGGTATGAAATCAACGGTTGAAAAACGTAATGCG
GCCATGTATTACCACAGCCTTTCCCATGTAATAACTGTGGGTGGCCTTTGCAGATGATACAGATCCT
CATGTCCATCACAGCACTCACCTAGGAGAAAAATCTTATAAATGTGACCAGTATGGAAGAAGCTTTAGT
CAGAGCCAAGATCTTATCGTTCATTGTAATAACTACTCTGGCAAGACTCCCTATGAATCCACGAATGG
CCTATGGGCTGCAAACAGAGCTCAGACCTTCCCAGATATCAGAAAGTCTCCTCAGGAGACAAACCCTAC
AAATGTAAGAATGTGGCAAGGGCTTCAGGCGCAGCTCTCCCTTCAACAACCATCATCGAGTCCACACA
GGGAGATGCCCTACAAATGCGATGAATGTGGAAAGGGTTTGGATTTAGGTCACCTCTTTGTATTTCAT
CAGGGAGTACACACAGGGAAAAAGCCCTATAAATGTGAAGAGTGTGGAAAGGGCTTTGATCAGAGCTCC
AACCTTCTGTCCATCAGAGAGTCCACACTGGAGAGAAGCCCTACAAATGCAGTGAGTGTGGCAAGTGC
TTAGTTCAAGCTCCGTTCTTCAAGTCCACTGGAGTTTACACAGGGGAGAAACCATATAGGTGTGGT
GAGTGTGGAAAGGGCTTCAGCCAATGTACACACCTTACATTACCAGAGAGTCCACACAGGGGAGAAA
CCATACAAATGCAATGTGTGGAAAGGATTTTGCATAGCTCTGTTCTTCCACTCATCAGAGAGTT
CACACTGGAGAAAAACCATATAAATGTGAAGTGTGGAAAGTGTTCAGTTACAGTTTCATATTTTAC
TTACATCAAAGAGATCACATCAGAGAGAAACCATATAAATGTGATGAGTGTGGTAAAGGCTTCAGCCGG
AATTCAGATCTTAATGTTACCTCAGAGTCCACACAGGAGAGAGGCCCTATAAGTGAAGGCATGTGGT
AAGGGCTTCAGTCGTAATTCATACCTCCTTGCCCATCAGAGAGTGCATATAGATGAGACACAGTACACA
CACTGTGAGCGTGGAAAGGACCTTCTGACTCATCAAAGACTACATGAGCAGAGAGAAACATTATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: Sgfl-Mlul

Plasmid Map:



ACCN: NM_001291490

Insert Size: 1308 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001291490.1](#)

RefSeq Size: 2841 bp

RefSeq ORF: 1308 bp

Locus ID: 26974

UniProt ID: [Q96NJ3](#)

Cytogenetics: 19q13.31

Protein Families: Transcription Factors

MW: 50 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (4) uses an alternate splice site in the 5' region and initiates translation at a downstream in-frame start codon, compared to variant 1. The encoded isoform (c) has a shorter N-terminus than isoform a. Variants 4 and 5 encode the same isoform (c). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.