

Product datasheet for SC336126

IFIT3 (NM_001289758) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	IFIT3 (NM_001289758) Human Untagged Clone
Tag:	Tag Free
Symbol:	IFIT3
Synonyms:	CIG-49; cig41; GARG-49; IFI60; IFIT4; IRG2; ISG60; P60; RIG-G
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336126 representing NM_001289758. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTACAACCTGTTGGCCTACATAAAACACCTAGATGGTAAACGAGGCAGCCCTGGAATGCTTACGG
CAAGCTGAAGAGTTAATCCAGCAAGAACATGCTGACCAAGCAGAAATCAGAAGTCTAGTCACTTGGGGA
AACTACGCTGGGTCTACTATCACTTGGGCAGACTCTCAGATGCTCAGATTTATGTAGATAAGGTGAAA
CAAACCTGCAAGAAATTTCAAATCCATACAGTATTGAGTATTCTGAACCTGACTGTGAGGAAGGGTGG
ACACAACCTGAAGTGTGAAGAAATGAAAGGGCGAAGGTGTGTTTTGAGAAGGCTCTGGAAGAAAAGCCC
AACACCAGAAATTCCTCTGGACTGGCAATTGCGATGTACCATCTGGATAATCACCCAGAGAAACAG
TTCTCTACTGATTTTTGAAGCAGGCCATTGAGCTGAGTCTGATAACCAATACGTCAAGTTCTCTTG
GGCCTGAAACTGCAGAAGATGAATAAAGAGCTGAAGGAGAGCAGTTTGTGAAGAAGCCTTGAAAAG
TCTCCTTGCCAAACAGATGTCCTCCGAGTGCAGCCAAATTTTACAGAAGAAAAGGTGACCTAGACAAA
GCTATTGAACTGTTTCAACGGGTGTTGGAATCCACACCAAAACAATGGCTACCTCTATCACCAGATTGGG
TGCTGCTACAAGGCAAAAGTAAGACAAATGCAGAATACAGGAGAACTGAAGCTAGTGGAAATAAAGAG
ATGATTGAAGCACTAAAGCAATATGCTATGGACTATTCGAATAAAGCTCTTGAGAAGGGACTGAATCCT
CTGAATGCATACTCCGATCTCGTGAATTCCTGGAGACGGAATGTTATCAGACACCATTCAATAAGGAA
GTCCCTGATGCTGAAAAGCAACAATCCCATCAGCGCTACTGCAACCTTCAGAAATAAATGGGAAGTCT
GAAGACACTGCTGTGCAACATGGTTTAGAGGGTTTGTCCATAAGCAAAAAATCAACTGACAAGGAAGAG
ATCAAAGACCAACCACAGAATGTATCTGAAAATCTGCTTCCACAAAAATGCACCAAAATTATTGGTATCTT
CAAGGATTAATTCATAAGCAGAATGGAGATCTGCTGCAAGCAGCCAAATGTTATGAGAAGGAACTGGGC
CGCCTGCTAAGGGATGCCCTTCAGGCATAGGCAGTATTTTCTGTGAGCATCTGAGCTTGAGGATGGT
AGTGAGGAAATGGGCCAGGGCGCAGTCCAGTCCAGTCCAGAGAGCTCCTCTCTAACTCAGAGCAACTG
AACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001289758

Insert Size: 1317 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001289758.1](#)

RefSeq Size: 2696 bp

RefSeq ORF: 1317 bp

Locus ID: 3437

Cytogenetics: 10q23.31

MW: 49.8 kDa

Gene Summary:

IFN-induced antiviral protein which acts as an inhibitor of cellular as well as viral processes, cell migration, proliferation, signaling, and viral replication. Enhances MAVS-mediated host antiviral responses by serving as an adapter bridging TBK1 to MAVS which leads to the activation of TBK1 and phosphorylation of IRF3 and phosphorylated IRF3 translocates into nucleus to promote antiviral gene transcription. Exhibits an antiproliferative activity via the up-regulation of cell cycle negative regulators CDKN1A/p21 and CDKN1B/p27. Normally, CDKN1B/p27 turnover is regulated by COPS5, which binds CDKN1B/p27 in the nucleus and exports it to the cytoplasm for ubiquitin-dependent degradation. IFIT3 sequesters COPS5 in the cytoplasm, thereby increasing nuclear CDKN1B/p27 protein levels. Upregulates CDKN1A/p21 by downregulating MYC, a repressor of CDKN1A/p21. Can negatively regulate the apoptotic effects of IFIT2.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (3) contains an alternate exon compared to variant 1. The resulting isoform (b) uses a translation initiation site at a downstream AUG that results in an isoform (b) with a shorter N-terminus compared to isoform 1. Variants 3 and 4 both encode the same isoform (b). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.