

Product datasheet for **SC336153**

ACPL2 (PXYLP1) (NM_001282728) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ACPL2 (PXYLP1) (NM_001282728) Human Untagged Clone
Tag:	Tag Free
Symbol:	PXYLP1
Synonyms:	ACPL2; HEL124; XYLP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336153 representing NM_001282728. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTAGCAAGAGTCGAAAGAGAATCATGCCGACCCTGTGACGGAGCCCCCTGTGACAGACCCCGTT
TATGAAGCTCTTTTGTACTGCAACATCCCCAGCGTGGCCGAGCGCAGCATGGAAGGTCATGCCCCGCAT
CATTTTAAAGCTGGTCTCAGTGATGTGTTCAATCGCCACGGAGACAGGTACCCACTGTATGTCATTCCC
AAAACAAAGCGACCAGAAATTGACTGCACTCTGGTGGCTAACAGGAAACCGTATCACCCAAAACCTGGAA
GCTTTTATTAGTCACATGTCAAAGGATCCGGAGCCTCTTTCGAAAGCCCCTTGAACCTCTTGCCTCTT
TACCCAAATCACCAATTGTGTGAGATGGGAGAGCTCACACAGACAGGAGTTGTGCAGCATTGTCAGAAC
GGTCAGCTGCTGAGGGATATCTATCTAAAGAAACACAACTCCTGCCCAATGATTGGTCTGCAGACCAG
CTCTATTTAGAGACCACTGGGAAAAGCCGACCCTACAAAGTGGGCTGGCCTTGCTTTATGGCTTTCTC
CCAGATTTTACTGGAAGAAGATTTATTTAGGCACCAGCCAAGTGCCTGTTCTGCTCTGGAAGCTGC
TATTGCCCGGTAAGAAACAGTATCTGGAAAAGGAGCAGCGTCGTGAGTACCTCTACGTTTGAAAAAC
AGCCAGCTGGAGAAGACCTACGGGGAGATGGCCAAAGATCGTGGATGTCCCACCAAGCAGCTTAGAGCT
GCCAACCCATAGACTCCATGCTCTGCCACTTCTGCCACAATGTCAGCTTCCCTGTACCAGAAATGGC
TGTGTTGACATGGAGCACTTCAAGGTAATTAAGACCCATCAGATCGAGGATGAAAGGGAAGACGGGAG
AAGAAATTGACTTCGGGTATTCTCCTGGGTGCCACCCATCCTGAACCAAAACCATCGGCCGGATG
CAGCGTGCCACCGAGGGCAGAAAAGAAGAGCTCTTTGCCCTCTACTCTGCTCATGATGCACTCTGTCA
CCAGTTCTCAGTGCCCTTGGCCCTTTTCAAGGCCAGGTTCCCAAGGTTTGCAGCCAGGTTGATCTTTGAG
CTTTGGCAAGACAGAGAAAAGCCAGTGAACATTCCTGTCGGATTCTTTACAATGGCGTCGATGTACA
TTCCACACCTCTTTCTGCCAAGACCACCACAAGCGTTCTCCCAAGCCATGTGCCCGCTTGAAAACCTTG
GTCCGCTTTGTGAAAAGGACATGTTTGTAGCCCTGGGTGGCAGTGGTACAAATTATTATGATGCATGT
CACAGGAAGGATTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001282728

Insert Size: 1329 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282728.1](#)

RefSeq Size: 3379 bp

RefSeq ORF: 1329 bp

Locus ID: 92370

UniProt ID: [Q8TE99](#)

Cytogenetics: 3q23

Protein Families: Transmembrane

MW: 51 kDa

Gene Summary: Responsible for the 2-O-dephosphorylation of xylose in the glycosaminoglycan-protein linkage region of proteoglycans thereby regulating the amount of mature glycosaminoglycan (GAG) chains. Sulfated glycosaminoglycans (GAGs), including heparan sulfate and chondroitin sulfate, are synthesized on the so-called common GAG-protein linkage region (GlcUA β 1-3Gal β 1-3Gal β 1-4Xyl β 1-O-Ser) of core proteins, which is formed by the stepwise addition of monosaccharide residues by the respective specific glycosyltransferases. Xylose 2-O-dephosphorylation during completion of linkage region formation is a prerequisite for the initiation and efficient elongation of the repeating disaccharide region of GAG chains.
[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (3) differs in its 5' UTR and initiates translation at a downstream start codon, compared to variant 1. The encoded isoform (2) has a shorter N-terminus, compared to isoform 1.