

Product datasheet for SC336176

ANKRD2 (NM_001291218) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ANKRD2 (NM_001291218) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANKRD2
Synonyms:	ARPP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336176 representing NM_001291218. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC CGCATCGCC
ATGCAGCCAGCAGTTCCCTTTTGC AACTCAGTTGCCTTGGCTAAGGCCCAACCAGTCCCAAGTGCTCAG
GACAGAGGCTGGACGGGCTCCTTCTGCACCCCTGCTCTTGGGACACAGTGCCCTCCGGCTCTAATAGG
CCAGGAGTTGGGGGGCAACTGGCTCTGCTCCCTGGCCCTGGCTCCCCCTGCTCCCTGGGGCTGCCCGA
GGTGAAGGTGACAGGTGGGGGAGGCAGGTGGAGAATTGGCCAGTGAGCTCATGGCAAAGGCGCCAGC
TGGGCAGGGTGGTCTCTGGCCTATAAAGCCCCGAGGCCCTGTGGCCTGCAGAGCGGTTATGGAC
GGCACCATGGAGGACTCCGAGGCGGTGCAGAGGGCCACAGCGCTCATCGAGCAGCGGCTGGCACAGGAG
GAGGAGAATGAGAACTCCGAGGAGACGCACGCCAGAAGCTGCCCATGGACTTGTGGTCTGGAGGAT
GAGAAGCACCCAGGGGCTCAGAGTGCAGCCCTGCAGAAGGTGAAGGGCCAAGAGCGCGTGCGAAGACG
TCCTTGGACCTGCGGCGGGAGATCATCGATGTGGGCGGGATCCAGAACCTCATCGAGCTGCGGAAGAAA
CGCAAGCAGAAGAAGCGGGACGCTCTGGCCGCCTCGCATGAGCCGCCCCAGAGCCCAGGAGATCACT
GGCCCTGTGGATGAGGAGACCTTCTGAAAGCTGCGGTGGAGGGGAAAATGAAGGTCATTGAGAAGTTC
CTGGCTGACGGGGGTCAGCCGACACGTGCGACCAGTTCCTGCGACAGCACTGCACCGAGCTTCCCTG
GAAGGCCACATGGAATCCTGGAGAAGCTTCTAGATAATGGGGCCACTGTGGACTCCAGGATCGGCTG
GACTGCACAGCCATGCATTGGGCTGCCGCGGGGCCACTTAGAGGTGGTAAAATCTGCAAAGCCAT
GGAGCAGACACCAATGTGAGGGATAAGCTGCTGAGCACCCGCTGCACGTGGCAGTCCGGACAGGGCAG
GTGGAGATTGTGGAGCACTTTCTATCCCTGGGCTGGAATCAATGCCAGAGACAGGGAAGGGGATACT
GCCCTGCATGACGCTGTGAGGCTCAACCCTACAAAATCATCAAAGTGTCTCTGCATGGGGCTGAC
ATGATGACCAAGAACCTGGCAGGAAAGACCCCGACGGACCTGGTGCAGCTCTGGCAGGCTGATACCCGG
CACGCCCTGGAGCATCCTGAGCCGGGGCTGAGCATAACGGGCTGGAGGGGCTAATGATAGTGGGCGA
GAGACCCCTCAGCTGTGCCAGCCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
```



Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001291218

Insert Size: 1341 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001291218.1](#)

RefSeq Size: 1520 bp

RefSeq ORF: 1341 bp

Locus ID: 26287

UniProt ID: [Q9GZV1](#)

Cytogenetics: 10q24.2

MW: 48.9 kDa

Gene Summary:

This gene encodes a protein that belongs to the muscle ankyrin repeat protein (MARP) family. A similar gene in rodents is a component of a muscle stress response pathway and plays a role in the stretch-response associated with slow muscle function. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Mar 2014]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longest isoform (c, also known as L-Ankrd2).