

## Product datasheet for **SC336178**

### LILRB4 (NM\_001278427) Human Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	LILRB4 (NM_001278427) Human Untagged Clone
Tag:	Tag Free
Symbol:	LILRB4
Synonyms:	CD85K; ILT-3; ILT3; LIR-5; LIR5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336178 representing NM_001278427. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATCCCCACCTTCACGGCTCTGCTCTGCCTCGGGCTGAGTCTGGGCCCCAGGACCCACATGCAGGCA
GGGCCCTCCCCAAACCCACCCTCTGGGCTGAGCCAGGCTCTGTGATCAGCTGGGGAACTCTGTGACC
ATCTGGTGTGAGGGACCTGGAGGCTCGGGAGTACCGTCTGGATAAAGAGGAAAGCCAGCACCCCTGG
GACAGACAGAACCCTGGAGCCCAAGAACAAGGCCAGATTCTCCATCCCATCCATGACAGAGGACTAT
GCAGGGAGATACCGCTGTTACTATCGCAGCCCTGTAGGCTGGTCACAGCCAGTGACCCCTGGAGCTG
GTGATGACAGGAGCCTACAGTAAACCCACCCTTTCAGCCCTGCCGAGTCTCTTGTGACCTCAGGAAAG
AGCGTGACCCTGCTGTGTGATCAGCCAGCCCAATGGACACTTTTCTTCTGATCAAGGAGCGGGCAGCC
CATCCCCACTGTCATCTGAGATCAGAGCACGGAGCTCAGCAGCACCCAGGCTGAATCCCCATGAGTCTT
GTGACCTCAGTGCACGGGGGACCTACAGGTGCTTTCAGCTCACACGGCTTCTCCACTACCTGCTGTCA
CACCCAGTGACCCCTGGAGCTCATAGTCTCAGGATCCTTGGAGGGTCCCAGGCCCTCACCCACAAGG
TCCGTCTCAACAGCTGCAGGCCCTGAGGACCAGCCCTCATGCCTACAGGGTCACTCCCCACAGTGGT
CTGAGAAGGCACTGGGAGTACTGATCGGGTCTTGGTGGTCTCCATCCTGCTTCTCCCTCCTCCTC
TTCTCCTCCTCAACACTGGCGTCAGGAAAACACAGGACATTGGCCCAGAGACAGGCTGATTTCAA
CGTCTCCAGGGCTGCCGAGCCAGAGCCCAAGGACGGGGCTACAGAGGAGGTCAGCCAGCCAGCTGCT
GACGCTCAGGGGAAAATCTGTGCTGCCGTGAAGAACACACAGCCTGAGGACGGGGTGGAAATGGAC
ACTCGGAGCCACACGATGAAGACCCCCAGGCAGTGACGTATGCCAAGGTGAAACTCCAGACCTAGG
AGAGAAATGGCCTCTCCTCCCTCCCACTGTCTGGGGAATTCCTGGACACAAAGGACAGACAGGCAGAA
GAGGACAGACAGATGGACTGAGGCTGCTGCATCTGAAGCCCCCAGGATGTGACCTACGCCCGGCTG
CACAGCTTTACCCTCAGACAGAAGGCAACTGAGCCTCCTCCATCCCAGGAAGGGGCTCTCCAGCTGAG
CCCAGTGTCTATGCCACTCTGGCCATCCACTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
```



Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM\_001278427

Insert Size: 1344 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001278427.3](#)

RefSeq Size: 3899 bp

RefSeq ORF: 1344 bp

Locus ID: 11006

UniProt ID: [Q8NHJ6](#)

Cytogenetics: 19q13.42

Protein Families: Transmembrane

**MW:** 49.2 kDa

**Gene Summary:** This gene is a member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at chromosomal region 19q13.4. The encoded protein belongs to the subfamily B class of LIR receptors which contain two or four extracellular immunoglobulin domains, a transmembrane domain, and two to four cytoplasmic immunoreceptor tyrosine-based inhibitory motifs (ITIMs). The receptor is expressed on immune cells where it binds to MHC class I molecules on antigen-presenting cells and transduces a negative signal that inhibits stimulation of an immune response. The receptor can also function in antigen capture and presentation. It is thought to control inflammatory responses and cytotoxicity to help focus the immune response and limit autoreactivity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]  
Transcript Variant: This variant (2) uses an alternate in-frame splice site in the 3' coding region, compared to variant 1. The encoded isoform (2) is shorter, compared to isoform 1.  
Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.