

Product datasheet for **SC336181**

TTC23 (NM_001288616) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TTC23 (NM_001288616) Human Untagged Clone
Tag:	Tag Free
Symbol:	TTC23
Synonyms:	HCC-8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336181 representing NM_001288616. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAAGAATCACAGGAAACCCACATATCCAACCACCTAGATGAAGTTGTTGCTGCTGTTAGCATCACT
CATAGAAAAGAAGTTCCAAAACAAGCTGCTTCAGACAGCACTATCCAGCCTCCTCGAGAGAACTCCAC
CTCTGTGAAGAGAAAGCAAAGTCTATTCCAACAGTCATGAGTACAAACAGGCCGTCATGAGCTTGTG
CGTTGCGTAGCACTGACAAGAATTTGCTATGGAGACTCACATTGGAACTAGCAGAGGCACATGTTAAT
CTGGCTCAAGGCTACCTCCAGCTGAAAGGACTGTCAGTCAAGCAAAACAACATGCAGAAAAAGCCAGA
CAAATCCTCGCCAACTCCATTGTGCTCCCTATAGTGAGAATACAGATGTTTTCAAGTTTTCCATTGAG
CTTTCCATACCATGGGCAGAGCTTTACTCTCCCTTCAAAAATTTAAGGAAGCTGCAGAGAAATTTGACA
AAAGCAGAGAGACTTTCAAAGGAGCTGCTACAATGTGGAAGAATTATAAAGGAAGAATGGATAGAAATT
GAAGCACGGATCAGATTATCATTGACACAGGTGTATCAAGGTGAGAAAGTCAAAGAAGCTTTGTCC
CACTATCAAGCAGCTTTGGAATATGTTGAGATCAGTAAAGGTGAAACAAGTCGTGAGTGTGTACCCATA
TTGAGAGAATTAGCAGGTGTAGAGCAAGCCCTGGGACTCCACGATGTATCCATCAACCCTTCTCCAG
GCACATCTTATCATCCTGAGTAGAAGCCCTCTCAAGTGGAGGCAGCAGACTCGGCACACATCGTCGCC
CATGCTGCTGTCGCTTCAGGGAGACACGAGCACCATGATGTAGCTGAGCAGTATTTTCAAGAGAGCATG
GCTCATCTTAAGGATTCTGAAGGGATGGGAAGAACCAATTTCTTTCAATTCAAGATGAATTTGGCAT
TTTCTACAAATGACTGGACAAAAAGAGAGAGCAACCTCGATCCTGAGAGAGTCCCTGGAAGCCAAAGTG
GAAGCATTGGCGATTTAGTCCCGAGGTGGCAGAGACATACCGCTCCTGGGAGGAGCAGACCTGGCG
CAGGGGAACCACAGTGGGCCCCGAAGAACTGAAGAAGTGTCTCCAGATCCAGACCCTCTTATATGGA
CCGCAGGACAAAAGGACTCTGGCCACCCAGCAGGCCATGGGCATGCTGTCCACGGCCCCAAGGTTGCT
TCGAAGCCAAGGCAGGCATCAAAGCCAAAGTGGCCTTCTGCACCAGCATCCCTCAGGACACCCTGCTG
GGGAAGGCCCGCCCGGCACAACAGCAGACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001288616

Insert Size: 1344 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001288616.1](#)

RefSeq Size: 3324 bp

RefSeq ORF: 1344 bp

Locus ID: 64927

UniProt ID: [Q5W5X9](#)

Cytogenetics: 15q26.3

MW: 50 kDa

Gene Summary:

Participates positively in the ciliary Hedgehog (Hh) signaling.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (9) differs in the 5' UTR, compared to variant 8. Variants 8 and 9 encode the same protein.