

Product datasheet for **SC336186**

LILRB4 (NM_001278426) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LILRB4 (NM_001278426) Human Untagged Clone
Tag:	Tag Free
Symbol:	LILRB4
Synonyms:	CD85K; ILT-3; ILT3; LIR-5; LIR5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336186 representing NM_001278426. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATCCCCACCTTACGGCTCTGCTCTGCCTCGGGCTGAGTCTGGGCCCCAGGACCCACATGCAGGCA
GGGCCCTCCCCAAACCCACCCTCTGGGCTGAGCCAGGCTCTGTGATCAGCTGGGGAACTCTGTGACC
ATCTGGTGTGAGGGACCTGGAGGCTCGGGAGTACCGTCTGGATAAAGAGGAAAGCCAGCACCCCTGG
GACAGACAGAACCCTGGAGCCCAAGAACAAGGCCAGATTCTCCATCCCATCCATGACAGAGGACTAT
GCAGGGAGATACCGCTGTTACTATCGCAGCCCTGTAGGCTGGTCACAGCCAGTGACCCCTGGAGCTG
GTGATGACAGGAGCCTACAGTAAACCCACCCTTTCAGCCCTGCCGAGTCCCTTTGTGACCTCAGGAAAG
AGCGTGACCCTGCTGTGTGATCAGTCACGGAGCCCAATGGACTTTTTCTTCTGATCAAGGAGCGGGCAGCC
CATCCCCACTGTCATCTGAGATCAGAGCACGGAGCTCAGCAGCACCCAGGCTGAATTCCCCATGAGTCCT
GTGACCTCAGTGCACGGGGGACCTACAGGTGCTTTCAGCTCACACGGCTTCTCCACTACCTGCTGTCA
CACCCAGTGACCCCTGGAGCTCATAGTCTCAGGATCCTTGGAGGGTCCCAGGCCCTCACCCACAAGG
TCCGTCTCAACAGCTGCAGGCCCTGAGGACCAGCCCTCATGCCTACAGGGTCACTCCCCACAGTGGT
CTGAGAAGGCACTGGGAGTACTGATCGGGTCTTGGTGGTCTCCATCCTGCTTCTCCTCCTCCTCCTC
TTCTCCTCCTCAACACTGGCGTCAGGGAAAACACAGGACATTGGCCCAGAGACAGGCTGATTTCAA
CGTCTCCAGGGGCTGCCGAGCCAGAGCCCAAGGACGGGGCCTACAGAGGAGTCCAGCCAGCTGCT
GACGTCAGGGGAGAAAATTTCTGTGCTGCCGTGAAGAACACACAGCCTGAGGACGGGGTGGAAATGGAC
ACTCGGCAGAGCCACACGATGAAGACCCCGAGGAGTACGTATGCCAAGGTGAAAACACTCCAGACCT
AGGAGAGAAATGGCTCTCCTCCTCCCACTGTCTGGGAAATTCCTGGACACAAAGGACAGACAGGCA
GAAGAGGACAGACAGATGGACTGAGGCTGCTGCATCTGAAGCCCCCAGGATGTGACCTACGCCCGG
CTGCACAGCTTTACCCTCAGACAGAAGGCAACTGAGCCTCCTCCATCCCAGGAAGGGGCTCCTCAGCT
GAGCCAGTGTCTATGCCACTCTGGCCATCCACTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
```



Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001278426
Insert Size:	1347 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001278426.3</u>
RefSeq Size:	3902 bp
RefSeq ORF:	1347 bp
Locus ID:	11006
UniProt ID:	<u>Q8NHJ6</u>
Cytogenetics:	19q13.42
Protein Families:	Transmembrane
MW:	49.3 kDa
Gene Summary:	<p>This gene is a member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at chromosomal region 19q13.4. The encoded protein belongs to the subfamily B class of LIR receptors which contain two or four extracellular immunoglobulin domains, a transmembrane domain, and two to four cytoplasmic immunoreceptor tyrosine-based inhibitory motifs (ITIMs). The receptor is expressed on immune cells where it binds to MHC class I molecules on antigen-presenting cells and transduces a negative signal that inhibits stimulation of an immune response. The receptor can also function in antigen capture and presentation. It is thought to control inflammatory responses and cytotoxicity to help focus the immune response and limit autoreactivity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) encodes isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>