

Product datasheet for **SC336188**

KIAA1958 (NM_001287038) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIAA1958 (NM_001287038) Human Untagged Clone
Tag:	Tag Free
Symbol:	KIAA1958
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336188 representing NM_001287038. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**GCGATCGCC**
ATGGAGGATTGTCTTCATACCTCATCTGAGAATCTGTCCAAATTGGTCAGCTGGGCCCATAGCCATGGG
ACTATTTCGAGCCTCATTCCAAACCTGAAACACTTGCTTTCTGAAGGTTCCCATGGGAACCTGACAGCA
ATGTGGGGCTGTAGTGCTGGCCATGCTTATCACTGGCCACTAACAGCTACTTGCAGAGCTGGGTCCCAA
GAGAGGGTCTGTTCCAGGATAACAGAAGTTTAACTCTGATAGTCCAGTATAATCGGGGTGCCCTCT
GAGACACAGACTAGCCCTGTTGAAAGGTACCCTGGGAGACCAGTGAAAGCAAAGCTAGACTGTAACCGG
ACCAGAGACTCTGTGACTTCTCCTACTGTAGTGAGCCCTCTGAACTGGATGAACTGTTGAAGAATAT
GAAGATGAGAACCCTGTTTGACATGGTTGTGAGTCTTCTGTTACAGATGAGGATAGTGACTTTGAA
CCCCAAACCCAAAGGCCTCAAAGCATTGCTCGCAAAAGACCTGGGGTAGTCCCATCTTCCCTCCATTCA
AGCTCCCAGACGCAGATGGTTGACGAATGCAGCAATGATGTCATCATCAAGAAAATCAAACAAGAAATC
CCCGAAGATTATTACATTGTGGCAATGCAGAACTGACAGGAGGAGTAGATGGACCAGCCCTGTCCTTG
ACACAGATGGCAAAACCAAGCCTCAGACTCACGCTGGTCCCTCCTGTGTAGGGTCTGCTAAACTGATT
CCCCATGTCACATCTGCCATCAGCACGGAGCTAGACCCACACGGTATGTCTGCATCCCCCTCTGTGATC
TCCAGACCAATTGTCCAGAAGACTGCTAGGGTATCTCTGGCTTACCAAACAGAGGACCCCTGGTACA
CATGGCACCAACCAACAGGTGGCCATGCAAATGCCTGTGAGCACATCCCATCCTAACAAACAGATCAGT
ATCCCCCTGTCTGCCCTGCAGCTGCCTGGACAGGATGAGCAAGTTGCCTCTGAAGAGTTCCTGTCCCAT
CTGCCCAGCCAGGTCTCCTCCTGTGAGGTAGCCCTTTCTCCCTCAGTTAACACAGAGCCAGAAGTGAGC
TCCAGTCAGCAGCAGCCCCAGTCGCTCCAGCCATAACCACTGAGGCCACAGCACAGTGCATACCAGCT
TATTCCACTAAGCTCAACAAATTTCTGTATTTAATATTAATGATGACTTGAATGATCTGTGTACCAGT
GCAGTAAGCCCAATACTACAAAGCCACGCGGTACGCTTGAATGTGTGGCGTTATTGGTGCATGACC
AACGGGCTCAAAGACCACACAGACATACCAAG**TGA**
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC

Restriction Sites: SgfI-MluI



View online »

ACCN:	NM_001287038
Insert Size:	1347 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001287038.1</u>
RefSeq Size:	7565 bp
RefSeq ORF:	1347 bp
Locus ID:	158405
UniProt ID:	<u>Q8N8K9</u>
Cytogenetics:	9q32
MW:	48.6 kDa