

Product datasheet for SC336206

Chimaerin 2 (CHN2) (NM_001293072) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	Chimaerin 2 (CHN2) (NM_001293072) Human Untagged Clone
Tag:	Tag Free
Symbol:	CHN2
Synonyms:	ARHGAP3; BCH; CHN2-3; RHOGAP3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336206 representing NM_001293072. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATGCTGAAGAATACCAGCCTCCTATATGGAAATCATACTTATATCAGTTACAGCAAGAGGCACCT
CGTCCCAAGAGAATCATTGTCTCGGGAGGTGGAAAACAGACAAAATATTATGGAAGAGAGTTTCAT
GGGATCATCTCTCGGGAGCAGGCGGATGAGCTTCTGGAGGCGTGGAGGGTGCCTACATCCTTAGAGAA
AGCCAGCGGCAACCAGGATGCTACACGCTGGCTCTCAGGTTTGGAAACCAGACCTTAACTACAGGCTC
TTCCACGACGGGAAACACTTTGTGGGTGAGAAGAGGTTTGAAGTCGATTTCATGATCTGGTGACAGATGGC
TTGATAACTGTACATAGAAACAAAAGCTGCCGAGTACATTTCAAAAATGACAACCTAACCCCATCTAT
GAACACATTGGATATGCCACCCTACTCAGAGAAAAGTATCCAGAAGGCTGAGCAGGTCTAAAAATGAA
CCAAGAAAAACAAACGTCACACATGAAGAACACACAGCGGTGGAAAAGATCTCCTCCCTGGTTCCGAAGG
GCTGCCCTCACACACAACGACAACCACTTCAATTATGAGAAGACACAACTTTAAGGTCCACACGTTTC
CGAGGCCACACTGGTGTGAATATTGTGCCAATTTTCATGTGGGGGCTCATCGCCCAAGGGGTCCGGTGC
TCAGACTGTGGATTGAACGTACACAAACAGTGTCCAAGCACGTTCCCAATGACTGCCAACCTGATCTC
AAGAGGATCAAGAAAGTACTGTTGTGACCTCACAACACTTGTGAAGGCTCACAACACTCAGAGACCC
ATGGTGGTAGACATATGCATTGGGAAATTAAGCAAGAGGATTAATAATCGGAAGGCCTTTACAGAGTC
TCTGGTTCACTGAACACATTGAAGATGTCAAAATGGCATTGACAGAGATGGTGAAAAGGCCGATATA
TCTGCCAATGTCTATCCAGACATAAACATCATCACTGGAGCCCTTAACTGTATTTTCAGAGACTTACCC
ATCCCTGTGCATACATATGATACCTATTCCAAATTTATAGATGCAGCAAAAATCTCCAATGCAGATGAG
AGGCTGGAAGCCGTCCATGAAGTGTGATGCTGCTGCCTCCTGCCACTATGAAACCCTCCGGTACCTA
ATGATCCACCTCAAAAAGGTTACTATGAATGAAAAGACAATTTTCATGAATGCAGAAAATCTGGGGATC
GTGTTTGGGCCCACTCTGATGAGGCCCTGAGGACAGCACCTGACCACCCTGCATGATATGCGGTAC
CAAAAGCTGATTGTGCAGATTTAATAGAAAACGAAGACGTTTTATTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001293072

Insert Size: 1362 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001293072.1](#)

RefSeq Size: 3218 bp

RefSeq ORF: 1362 bp

Locus ID: 1124

Cytogenetics: 7p14.3

MW: 52.6 kDa

Gene Summary:

This gene encodes a guanosine triphosphate (GTP)-metabolizing protein that contains a phorbol-ester/diacylglycerol (DAG)-type zinc finger, a Rho-GAP domain, and an SH2 domain. The encoded protein translocates from the cytosol to the Golgi apparatus membrane upon binding by diacylglycerol (DAG). Activity of this protein is important in cell proliferation and migration, and expression changes in this gene have been detected in cancers. A mutation in this gene has also been associated with schizophrenia in men. Alternative transcript splicing and the use of alternative promoters results in multiple transcript variants. [provided by RefSeq, May 2014]

Transcript Variant: This variant (6) differs in the 5' UTR and contains multiple additional 5' coding exons, compared to variant 1. It represents use of an alternate promoter and initiates translation from an alternate start codon. The encoded isoform (6) has a distinct N-terminus and is longer than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.