

Product datasheet for SC336222

Retinoic Acid Receptor beta (RARβ) (NM_001290216) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Retinoic Acid Receptor beta (RARβ) (NM_001290216) Human Untagged Clone
Tag:	Tag Free
Symbol:	RARB
Synonyms:	HAP; MCOPS12; NR1B2; RARbeta1; RRB2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336222 representing NM_001290216. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACCACCAGCGCCACGCATGTCGGTCCCAGCAGTGAACGGACACATGACTCACTATCCAGCCACA
CCCTACCCGTTACTCTTTCCACCTGTCTCGGAGGACTTCCCTGCCTCCCTCCATGGCCTTCATGGA
CATCCGCCTCCGAGTGGATGCAGCACCCCGTCGCCGGCAACAATTGAAACACAGAGCACCAGCTCTGAG
GAACTCGTCCCAAGCCCCCATCTCCACTTCTCCCCCTCGAGTGTACAAACCCTGCTTCGTCGCCAG
GACAAATCATCAGGGTACCACTATGGGGTCAGCGCCTGTGAGGGATGTAAGGGCTTTTTCCGCAGAAGT
ATTCAGAGAATATGATTTACACTTGTCCAGAGATAAGAAGTGTGTTATTAATAAAGTACCAGGAAT
CGATGCCAATACTGTGACTCCAGAAGTCTTGAAGTGGGAATGTCAAAGAATCTGTCAGGAATGAC
AGGAACAAGAAAAGAAGGAGACTTCGAAGCAAGAATGCACAGAGAGCTATGAAATGACAGCTGAGTTG
GACGATCTCACAGAGAAGATCCGAAAAGCTCACCAGGAACTTTCCCTTCACTCTGCCAGCTGGGTTAAA
TACACCACGAATTCAGTGTGACCATCGAGTCCGACTGGACCTGGGCCTCTGGGACAAATTCAGTGAA
CTGGCCACCAAGTGCATTATTAAGATCGTGGAGTTTGCTAAACGTCTGCCTGGTTTCACTGGCTTGACC
ATCGCAGACCAAAATACCTGCTGAAGGCCGCTGCCTGGACATCCTGATTCTTAGAATTTGACCAGG
TATACCCAGAACAAAGACACCATGACTTTCTCAGACGGCCTTACCCTAAATCGAAGTCAAGTGCACAA
GCTGGATTTGGTCTCTGACTGACCTTGTGTTCCCTTTGCCAACCAGCTCCTGCCTTTGGAAATGGAT
GACACAGAAACAGGCCTTCTCAGTGCCATCTGCTTAATCTGTGGAGACCGCCAGGACCTTGAGGAACCG
ACAAAAGTAGATAAGCTACAAGAACCATTGCTGGAAGCACTAAAAATTTATATCAGAAAAAGACGCC
AGCAAGCCTCACATGTTTCAAAGATCTTAATGAAATCACAGATCTCCGTAGCATCAGTGTAAAGGT
GCAGAGCGTGTAATTACCTTGAAGATGAAATTCCTGGATCAATGCCACCTCTCATTCAAGAAATGCTG
GAGAATTTCTGAAGGACATGAACCCTTGACCCCAAGTTCAAGTGGGAACACAGCAGAGCACAGTCTAGC
ATCTCACCAGCTCAGTGGAAAACAGTGGGGTCAGTCAGTACCCTCGTCAAATAA
ACGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001290216

Insert Size: 1368 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001290216.2](#)

RefSeq Size: 3265 bp

RefSeq ORF: 1368 bp

Locus ID: 5915

UniProt ID: [P10826](#)

Cytogenetics: 3p24.2

Protein Families: Druggable Genome, Nuclear Hormone Receptor, Transcription Factors

Protein Pathways:	Non-small cell lung cancer, Pathways in cancer, Small cell lung cancer
MW:	50.5 kDa
Gene Summary:	<p>This gene encodes retinoic acid receptor beta, a member of the thyroid-steroid hormone receptor superfamily of nuclear transcriptional regulators. This receptor localizes to the cytoplasm and to subnuclear compartments. It binds retinoic acid, the biologically active form of vitamin A which mediates cellular signalling in embryonic morphogenesis, cell growth and differentiation. It is thought that this protein limits growth of many cell types by regulating gene expression. The gene was first identified in a hepatocellular carcinoma where it flanks a hepatitis B virus integration site. Alternate promoter usage and differential splicing result in multiple transcript variants. [provided by RefSeq, Mar 2014]</p> <p>Transcript Variant: This variant (3, also known as beta-1) initiates from an upstream promoter and differs in the 5' UTR and 5' coding region compared to variant 1. The resulting isoform (3) has a longer and different N-terminus compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>