

Product datasheet for **SC336248**

MAP3K19 (NM_001282883) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAP3K19 (NM_001282883) Human Untagged Clone
Tag:	Tag Free
Symbol:	MAP3K19
Synonyms:	RCK; YSK4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336248 representing NM_001282883.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTTCTATGCCAAAACAGAAAGACATGCTGAGTCATTGCTTGACATTTGTCATGATACAAACTCT
TCTCCAAGTATTTGATGACAGTTACCAAAAATCAAAACATCATCTTGCAAAGCATCAGCAGAAGTGAG
GAGTTTCGACCAAGATGGTGACTGCAGTCAATCCACACTGGTAAATGAAGAAGAAGATCCCAGTGGTGGT
AGACAGGACTGGCAACCCAGGACAGAAGGTGTTGAGATCACTGTAACCTTTTCCAAGAGATGTCAGTCCT
CCCCAAGAAATGAGCCAAGAAGACTTAAAAGAAAAGAATCTGATAAACTCATCGTTCAAGAATGGGCA
CAAGCACATGCAGTTTCTCATCCAAATGAAATAGAAACGGTGGAGCTCAGGAAAAAGAAGCTGACCATG
CGGCCCTTAGTTTTGCAAAAAGAGGAAAGTTCAGGGAGCTCTGCAATGTGAACCTGGGCTTTTTGCTA
CCAAGTCTTGTTAGAAGTGAACATTTCCAAGTCTGTAACCAGAGAAGATGCTCCTCATTTTCTGAAG
GAGCAGCAAAGAAAATCTGAAGAGTTTTCGACCTCTCATATGAAGTACAGTGGCCGAAGCATCAAGGTA
TACTGTGGTCTCACTAGTCAAGGACAGCTAATAGCTGTA AACAGGTGGCTTTGGATACCTCTAATAAA
TTAGCTGCTGAAAAGGAATACCGGAAACTACAGGAAGAAGTAGATTTGCTCAAAGCACTGAAACATGTC
AACATTGTGGCCTATTTGGGGACATGCTTGCAAGAGAACACTGTGAGCATTTTTCATGGAGTTTGTCTCT
GGTGGCTCAATCTCTAGTATTATAAACCGTTTTGGGCCATTGCCTGAGATGGTGTCTGTAATATACG
AAACAAATACTTCAAGGTGTTGCTTATCTCCATGAGAAGTGTGTGGTACATCGCGATATCAAAGGAAAT
AATGTTATGCTCATGCCAACTGGAATAATAAGCTGATTGACTTTGGCTGTGCCAGGCGTTTTGGCCTGG
GCAGGTTTAAATGGCACCCACAGTGACATGCTTAAGTCCATGCATGGGACTCCATATTGGATGGCCCA
GAAGTCATCAATGAGTCTGGCTATGGACGGAATCAGATATCTGGAGCATTGGTTGACTGTGTTGAG
ATGGCTACAGGGAAGCCTCCACTGGCTTCCATGGACAGGATGGCCGCCATGTTTTACATCGGAGCACAC
CGAGGGCTGATGCCTCCTTTACCAGACCACTTCTCAGAAAATGCAGCAGACTTTGTGCGCATGTGCCTG
ACCAGGGACCAGCATGAGCGACCTTCTGCTCTCCAGCTCCTGAAGCACTCCTTCTTGAGAGAAGTAC
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001282883

Insert Size: 1383 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282883.1
RefSeq Size:	1783 bp
RefSeq ORF:	1383 bp
Locus ID:	80122
UniProt ID:	Q56UN5
Cytogenetics:	2q21.3
Protein Families:	Druggable Genome, Protein Kinase
MW:	51.9 kDa