

## Product datasheet for **SC336250**

### DCP1A (NM\_001290206) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	DCP1A (NM_001290206) Human Untagged Clone
Tag:	Tag Free
Symbol:	DCP1A
Synonyms:	HSA275986; Nbla00360; SMAD4IP1; SMIF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC336250 representing NM\_001290206.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTGATGTGGTAGAAGAGGAGACACGGCGATCCCAGCAAGCTGCTCGGGACAAACAGAGTCCCAGC
CAGGCCAATGGCTGCAGCGACCACAGGCCATCGACATCCTGGAGATGCTGAGCAGAGCCAAGGATGAG
TATGAGAGGAATCAGATGGGTGACTCAAATATCTCCAGCCCTGGGTTACAGCCAAGCACTCAGCTCTCC
AATCTGGGAAGCACCGAGACTCTAGAAGAAATGCCCTCCGGGTACAGGATAAGTCTGCTCCATCTGGA
CACAAGCATCTGACGGTAGAAGAGTTATTTGGAACCTCTTTGCCAAAGGAACAACCAGCAGTTGTGGGT
CTGGATTGAGAAGAAATGGAGAGGTTGCCAGGAGATGCCTCCCAGAAAGAGCCCAATTCATTCTACCA
TTTCCCTTTGAGCAGTTAGGAGGAGCCCTCAATCAGAAACCCTGGGTGTCCCTTCTGCTGCCACCAT
TCAGTCCAGCCTGAAATCACCACCCGGTCTAATCACTCCAGCCTCCATCACACAGTCCAATGAAAAG
CATGCTCCAACCTACACAATCCCGTTGAGCCCTGTTCTCAGTCCCCTCTGCCAGCTGAAGCTCCTACT
GCACAGGTTCCCCCAGCTTACCTCGAAACAGCACCATGATGCAGGCAGTGAAGACCAGCCTAGACAG
AGGTCTCCACTCCTGAACCAGCCAGTCCCTGAGCTAAGCCATGCCAGTCTGATTGCCAACAGAGCCCC
TTCAGGGCCCCATTGAACGTGACGAACACAGCTGGCACATCCCTCCCAAGCGTTGATCTTCTCCAGAAA
CTCAGGTTGACCCACAGCATGACCAAATACAGACACAACCCTTGGGAAAGGTGCAATGGTAGCCAGC
TTTTCTCCGGCAGCTGGTCAGCTAGCCACACCTGAGAGCTTATAGAGCCTCCCTCTAAGACAGCAGCA
GCAAGAGTGGCGGCCCTCAGCCTCCCTGAGCAACATGGTGCTTGTCCCTTCACTATGCAGCAGAAC
CAGGATCCTGAAGTATTTGTGCAGCCTAAGGTGTTATCCAGTGCCATCCCGTTGCAGGGCCCCACTG
GTTACTGCAACGACCCTGCAGTGTCTTCACTGCTGGCCCAAGTGTTCAGCAGACAGTACA
AGATCTTCGGACCTTGAGAGGAAAGCCAGCTCCCCTTCTCCTCTAACTATTGGAACGCCAGAAAGTACAG
AGAAAGCCTTCCATTATTCTCAGCAAGTCTCAGCTCCAGGATACATTAATACATTAATAAAGAATGAT
TCCAGTTCCTCAGTACTTTCATGAAGTCTACTTGCAGGTTCTGACCAAGAACAAGACAACCACAAC
CTATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** Sgfl-MluI

**Plasmid Map:**



**ACCN:** NM\_001290206

**Insert Size:** 1386 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001290206.1</a>
<b>RefSeq Size:</b>	5899 bp
<b>RefSeq ORF:</b>	1386 bp
<b>Locus ID:</b>	55802
<b>UniProt ID:</b>	<a href="#">Q9NPI6</a>
<b>Cytogenetics:</b>	3p21.1
<b>Protein Families:</b>	Transcription Factors
<b>Protein Pathways:</b>	RNA degradation
<b>MW:</b>	49.3 kDa
<b>Gene Summary:</b>	<p>Decapping is a key step in general and regulated mRNA decay. The protein encoded by this gene is a decapping enzyme. This protein and another decapping enzyme form a decapping complex, which interacts with the nonsense-mediated decay factor hUpf1 and may be recruited to mRNAs containing premature termination codons. This protein also participates in the TGF-beta signaling pathway. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Feb 2014]</p> <p>Transcript Variant: This variant (4) lacks an alternate internal exon in the 5' region, and it thus differs in its 5' UTR and initiates translation at a downstream in-frame start codon, compared to variant 1. The encoded isoform (d) is shorter at the N-terminus, compared to isoform a.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>