

Product datasheet for **SC336257**

TBCCD1 (NM_001286749) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	TBCCD1 (NM_001286749) Human Untagged Clone
Tag:	Tag Free
Symbol:	TBCCD1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336257 representing NM_001286749. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACAGACCGACTTAAGTTCGGATTCCGGCTTCCACCGCTTCTTAGTGCTTTCAGTGGACACGCTACAG
TTTCTGCTCTTCTATACATTCAACAGTTGAACAAGGTCTCCCTAAGGACATCTTTGATTGGCGAAGAG
TGGCCAGTCCCAGAAACAAATCTCAGTCTCCTGACCTGACTGAAAAATCTAATTGTCATAATAAGAAC
TGGAAATGATTACAGTACCAAGCTTTTGTCTATGATCATCTGTCTGATCTCCTCGAGCTGCTTTTAGAT
CCAAAACAACACTCACTGCATCATTTCAATCAACCCATAGTAGTCTAGTGTCTCGAGAAGCTGTTGTGGCG
CTCAGTTCCTTATTGAAGGTACAATAAGTAGAGCCAGGAAGATCTATCCACTTCATGAACCTGCACTG
TGGCAACCACTGCATGCAGATAGTGGCTTCTCAAAGATCTCTAAGACTTTCTCTTTCTACAAAACGGAA
ACCTGGTTGAGGTCCTGTTTACTGAGGAAATCCATTTGGTACATCAGCTTGCCCTCAAGTCTGAAAAGAAA
TTGGCTTGGGCTCATCAAGTTGAAGGGACCACAAAAGAGCTAAGATTGCTTGTAACTCATGTGGCC
CCTAGGATGCACCGACTGGTAGTGATGAGCCAGGTTTACAAGCAGACACTGGCTAAGAGCTCAGACACT
CTGGCGGGGGCACATGTAAGATTATCGTTGCAACGAATCTTTTATATATCTGCTCTCTCCCTACGA
TCTGTGACAATTGAGAAGTGCAGGAATAGCATCTTTGTCTTGGCCCTGTAGGGACTACACTCACCTC
CACAGTTGTGACAATGTTAAAGTCATTGCTGTTGCCATCGTTTGTCCATCTCTTCTACAACAGGTTGC
ATCTTTCAGTTCCTACGCCTACAGCCCACTTATTCTCTCTGGGAACAGACAGTAACCTTTGCCCT
TTTCATACACATTACCAATGCTAGAGGACCATATGGCCAGGACTGGCCTTGCTACAGTGCTAAGTAT
TGGGATAATCCAATGGTTGTGTGCAGAGAGAACAGCGACACAAGAGTCTTCCAGCTTTTACCACCTTGT
GAATTCTATGATTTATTATCCCTTTGAAATGGAAGGGGACACAACAGAGATACCCGGGGGCTTCCA
TCTGTATATCAGAAAGCACTGGGTCAAAGAGAACAAGATACAGATCTGGCAGAAAACGTGAAGGAG
GCTCATTTGACAAAGGATCAAAGGAAGCAGTCCAGGTAAGTACTGGTAGAGAACAAGTTTTATGAATGGTTG
ATTAATACAGGACATCGCCAACAGCTGGACAGCCTTGTACCCCTGCAGCAGGCTCCAAACAAGCAGCT
GGATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_001286749

Insert Size: 1386 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286749.1](#)

RefSeq Size: 2701 bp

RefSeq ORF: 1386 bp

Locus ID: 55171

UniProt ID: [Q9NVR7](#)

Cytogenetics: 3q27.3

MW: 52.3 kDa

Gene Summary:

Plays a role in the regulation of centrosome and Golgi apparatus positioning, with consequences on cell shape and cell migration.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (3) lacks part of the 5' coding region, and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (b) has a shorter and distinct N-terminus, compared to isoform a.