

## Product datasheet for **SC336264**

### PPP2R5E (NM\_001282180) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PPP2R5E (NM_001282180) Human Untagged Clone
Tag:	Tag Free
Symbol:	PPP2R5E
Synonyms:	B56E; B56epsilon
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336264 representing NM\_001282180.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCCTCAGCACCAACTACTCCTCCATCAGTGGATAAAGTAGACGGATTTTCTCGGAAGTCCGTCAGA
AAAGCCAGACAGAAGAGGTCGCAAAGTTCTCACAGTTTAGTCTCAAGCAAGCCTATTGAGTTAACA
CCTCTGCCGCTGCTAAAAGACGTTCCATCCTCAGAGCAGCCTGAACTGTTCTAAAGAACTTCAGCAG
TGCTGTGTCATTTTTGACTTCATGGACACGCTATCTGATCTTAAAATGAAAGAATACAAGCGCTCCACT
CTTAATGAACTGGTGGACTACATTACAATAAGCAGAGGCTGTTTGACAGAGCAGACTTACCCTGAAGTA
GTTAGAATGGTATCTTGAATATATTCAGAACTCTCCCTCCTAGTGACAGCAATGAATTTGATCCAGAA
GAAGATGAACCTACCCTTGAGGCATCGTGGCCACACTTACAGCTTGATATGAATTTTTCATACGATTT
TTGAAAGCCAAGAATCCAACCCAGCATTGCCAAAAATATATAGATCAGAAATTTGTATTACAGCTT
CTGGAGCTATTTGACAGCGAAGACCCTCGGGAACGGGACTACTTAAAAACAGTCTTACACAGAATTTAT
GGCAAGTTTCTTGGTCTTAGAGCATTTATCCGAAAACAGATTAACAATATTTTTCTAAGGTTTGTAT
GAAACAGAACACTTCAATGGTGTAGCTGAACTGCTGGAATATTAGGAAGTATTATCAATGGCTTTGCT
TTACCTCTTAAGGCAGAACAACAACAGTTTCTGGTGAAAGTATTGATCCCTTTACACACTGTCAGGAGC
TTACTCTTCCATGCACAGCTGGCATTGTATAGTACAGTTTCTGGAGAAAGATCCTTCACTCACA
GAACCAGTTATTAGGGGGTTAATGAAATTTTGGCCTAAAACATGTAGTCAAAAAGAGGTCATGTTCCCT
GGGGAACCTGGAAGAAATATTGGATGTGATTGAACCTTCACAATTTGTTAAAATCCAAGAACCTTTGTTT
AAACAAATCGCCAAGTGTGTATCTAGCCCCATTTTTCAGGTGGCAGAAAGAGCACTCTATTATTGGAA
AATGAATACATCATGAGTTTGTAGAAGAAAACCTAACGTCATCCTTCCCATCATGTTTCCAGCCTT
TATAGGATTTCAAAGAACATTGGAATCCGGCTATTGTGGCGTTGGTGTACAATGTGTTGAAGGCATT
ATGAAAATGAACAGCACCATGTTTGGACGAGCTGACAGCCACATAACAATGAGAAAAAGAAAAGGAG
CGTGAAGAATTGTGAAAAAATTGGAGGATCTGGAGTTAAAGAGAGGTCTTAGACGTGATGGAATAATT
CCAACCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM\_001282180

Insert Size: 1389 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001282180.1</a>
<b>RefSeq Size:</b>	4108 bp
<b>RefSeq ORF:</b>	1389 bp
<b>Locus ID:</b>	5529
<b>UniProt ID:</b>	<a href="#">Q16537</a>
<b>Cytogenetics:</b>	14q23.2
<b>Protein Families:</b>	Druggable Genome, Phosphatase
<b>Protein Pathways:</b>	Oocyte meiosis, Wnt signaling pathway
<b>MW:</b>	54 kDa
<b>Gene Summary:</b>	<p>The protein encoded by this gene belongs to the phosphatase 2A regulatory subunit B family. Protein phosphatase 2A is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. The B regulatory subunit might modulate substrate selectivity and catalytic activity. This gene encodes an epsilon isoform of the regulatory subunit B56 subfamily. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Aug 2013]</p> <p>Transcript Variant: This variant (3) uses an alternate in-frame splice junction at an end of a 3' exon compared to variant 1. The resulting isoform (b) has the same N- and C-termini but is shorter compared to isoform a.</p>