

Product datasheet for **SC336267**

MAP3K19 (NM_001018047) Human Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | MAP3K19 (NM_001018047) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | MAP3K19 |
| Synonyms: | RCK; YSK4 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)


Fully Sequenced ORF: >SC336267 representing NM_001018047.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTTCTATGCCAAAACAGAAAGACATGCTGAGTCATTGCTTGACATTTGTCATGATACAAACTCT
TCTCCAAGTATTGATGACAGTTACCAAAAATCAAAACATCATCTTGCAAAGCATCAGCAGAAGTGAG
GAGTTCGACCAAGATGGTGACTGCAGTCAATCCACTGGTTAATGAAGAAGAAGATCCCAGTGGTGGT
AGACAGGACTGGCAACCCAGGACAGAAGGTGTTGAGATCACTGTAACCTTTTCCAAGAGATGTCAGTCCT
CCCCAAGAAATGAGCCAAGAAGACTTAAAAGAAAAGAATCTGATAAACTCATCGTTCAAGAATGGGCA
CAAGCACATGCAGTTTCTCATCCAAATGAAATAGAAACGGTGGAGCTCAGGAAAAAGAAGCTGACCATG
CGGCCCTTAGTTTTGCAAAAAGAGGAAAGTTCAGGGAGCTCTGCAATGTGAACCTGGGCTTTTTGCTA
CCAAGTCTTGTTAGAAGTGAACATTTCCAAGTCTGTAACCAGAGAAGATGCTCCTCATTTTCTGAAG
GAGCAGCAAAGAAAATCTGAAGAGTTTTCGACCTCTCATATGAAGTACAGTGGCCGAAGCATCAAGAGG
CATAGTAGTGGGCTCAGGATATATGACAGGGAGGAAAATTTCTCATCTCAAATGAAAAGAAGATATTT
TCTGAAAATAGTTTTAAAGTCTGAAGAACCTATCCTATGGACCAAGGTAGATTTGCTCAAAGCACTGAAA
CATGTCAACATTGTGGCCTATTTGGGGACATGCTTGAAGAGAACACTGTGAGCATTTTTCATGGAGTTT
GTTCTGGTGGCTCAATCTCTAGTATTATAAACCGTTTTGGGCCATTGCCTGAGATGGTGTCTGTAAA
TATACGAAACAAATACTCAAGGTGTGGCTTATCTCCATGAGAAGTGTGTGGTACATCGCGATATCAAA
GGAAATAATGTTATGCTCATGCCAACTGGAATAATAAAGCTGATTGACTTTGGCTGTGCCAGGCGTTTG
GCCTGGCAGGTTTAAATGGCACCCACAGTGACATGCTTAAGTCCATGCATGGGACTCCATATTGGATG
GCCCCAGAAGTCATCAATGAGTCTGGCTATGGACGGAAATCAGATATCTGGAGCATTGGTTGACTGTG
TTTGAGATGGCTACAGGGAAGCCTCCACTGGCTTCCATGGACAGGATGGCCGCCATGTTTTACATCGGA
GCACACCGAGGGCTGATGCCTCCTTTACCAGACCACTTCTCAGAAAATGCAGCAGACTTTGTGCGCATG
TGCTGACCAGGGACCAGCATGAGCGACCTTCTGCTCTCCAGCTCCTGAAGCACTCCTTCTTGAGAGA
AGTCACTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001018047

Insert Size: 1389 bp

| | |
|-------------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001018047.2</u> |
| RefSeq Size: | 1789 bp |
| RefSeq ORF: | 1389 bp |
| Locus ID: | 80122 |
| UniProt ID: | <u>Q56UN5</u> |
| Cytogenetics: | 2q21.3 |
| Protein Families: | Druggable Genome, Protein Kinase |
| MW: | 52.5 kDa |