

Product datasheet for **SC336278**

RHD (NM_001282870) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RHD (NM_001282870) Human Untagged Clone
Tag:	Tag Free
Symbol:	RHD
Synonyms:	CD240D; DIIIc; RH; Rh4; RH30; RHCED; RhDCw; RHDel; RHDVA(TT); RhII; RhK562-II; RhPI; RHPII; RHXIII
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336278 representing NM_001282870.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAGCTCTAAGTACCCGCGTCTGTCCGGCGCTGCCTGCCCTCTGGGCCCTAACACTGGAAGCAGCT
CTCATTCTCCTCTTCTATTTTTTTACCCACTATGACGCTTCCTTAGAGGATCAAAAGGGGCTCGTGGCA
TCCTATCAAGTTGGCCAAGATCTGACCGTGATGGCGGCCATTGGCTTGGGCTTCCTCACCTCGAGTTTC
CGGAGACACAGCTGGAGCAGTGTGGCCTTCAACCTCTTCATGCTGGCGCTTGGTGTGCAGTGGGCAATC
CTGCTGGACGGCTTCCTGAGCCAGTTCCTTCTGGGAAGGTGGTCATCACACTGTTCAAGTATTGGCTG
GCCACCATGAGTGCTTTGTGGTGCTGATCTCAGTGGATGCTGTCTTGGGGAAGGTCAACTGGCGCAG
TTGGTGGTGTGGTGTGGTGGAGGTGACAGCTTAGGCAACCTGAGGATGGTCATCAGTAATATCTTC
AACACAGACTACCACATGAACATGATGCACATCTACGTGTTCCGACGCTATTTGGGCTGTCTGTGGCC
TGGTGCCTGCCAAAGCCTTACCCGAGGGAACGGAGGATAAAGATCAGACAGCAACGATACCCAGTTTG
TCTGCCATGCTGGCGCCCTCTTCTGTGGATGTTCTGGCCAAGTTTCAACTCTGCTCTGCTGAGAAGT
CCAATCGAAAGGAAGAATGCCGTGTTCAACACCTACTATGCTGTAGCAGTCAGCGTGGTGACAGCCATC
TCAGGGTCATCCTTGGCTCACCCCAAGGGAAGATCAGCAAGACTTATGTGCACAGTGGCGTGTGGCA
GGAGGCGTGGCTGTGGGTACCTCGTGTCACCTGATCCCTTCTCCGTGGCTTGGCATGGTGTGGTCTT
GTGGCTGGGCTGATCTCCGTCCGGGGAGCCAAGTACCTGCCGGGGTGTGTAACCGAGTGTGGGGATT
CCCCACAGCTCCATCATGGGCTACAACCTCAGCTTGTGGGTCTGCTGGAGAGATCATCTACATTGTG
CTGCTGGTGTGATACCGTCCGAGCCGGCAATGGCATGATTGGCTTCCAGGTCCTCCTCAGCATTGGG
GAACTCAGCTGGCCATCGTATAGCTCTACGCTCTGGTCTCCTGACAGTTTCTCATTGGCTGTTGG
ATTTTAAGCAAAAGCATCCAAGAAAAACAAGCCTGTTCAAAAACAAGACAACCTCCTCACTGTTGC
CTGCATTTGTACGTGAGAAACGCTCATGACAGCAAGTCTCCAATGTTCCGCGCAGGCACTGGAGTCAGA
GAAATGGAGTTGAATCCTTCTCTGCCACTCTTTGAGGAGAATCTACCATTTATTATGCACTGTAGA
ATACAACAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN: NM_001282870

Insert Size: 1392 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282870.1
RefSeq Size:	2763 bp
RefSeq ORF:	1392 bp
Locus ID:	6007
Cytogenetics:	1p36.11
Protein Families:	Transmembrane
MW:	50.2 kDa

Gene Summary: The Rh blood group system is the second most clinically significant of the blood groups, second only to ABO. It is also the most polymorphic of the blood groups, with variations due to deletions, gene conversions, and missense mutations. The Rh blood group includes this gene, which encodes the RhD protein, and a second gene that encodes both the RhC and RhE antigens on a single polypeptide. The two genes, and a third unrelated gene, are found in a cluster on chromosome 1. The classification of Rh-positive and Rh-negative individuals is determined by the presence or absence of the highly immunogenic RhD protein on the surface of erythrocytes. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (6) lacks an alternate coding exon compared to variant 1, that causes a frameshift. The resulting isoform (6) has a longer and distinct C-terminus compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.