

Product datasheet for **SC336280**

beta V Tubulin (TUBB) (NM_001293212) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	beta V Tubulin (TUBB) (NM_001293212) Human Untagged Clone
Tag:	Tag Free
Symbol:	TUBB
Synonyms:	CDCBM6; CSCSC1; M40; OK/SW-cl.56; TUBB1; TUBB5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336280 representing NM_001293212.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGAGATGTGGGGCCGAGGTGGGCGAATTGGGATCCCTCCAGGTCAGGGGTTTCGAGACCATCCTGGG
CAACAAAGCGAGACCCTCCCCATGCCAGTTTCTACAAAAAATAAAATTCTGGGAGGTGATCAGTGAT
GAACATGGCATCGACCCACCGGCACCTACCACGGGGACAGCGACCTGCAGCTGGACCGCATCTCTGTG
TACTACAATGAAGCCACAGGTGGCAAATATGTTCTCGTGCCATCCTGGTGGATCTAGAACCTGGGACC
ATGGACTCTGTTTCGCTCAGGTCCTTTTGGCCAGATCTTTAGACCAGACAACCTTTGATTTGGTCAGTCT
GGGGCAGGTAACAACGGCCAAAGGCCACTACACAGAGGGCGCCGAGCTGGTTGATTCTGTCTGGAT
GTGGTACGGAAGGAGGCAGAGAGCTGTGACTGCCTGCAGGGCTTCCAGCTGACCCACTACTGGGCGGG
GGCACAGGCTCTGAATGGGCACTCTCCTTATCAGCAAGATCCGAGAAGAATACCCTGATCGCATCATG
AATACCTTCAGTGTGGTGCCTTACCCAAAGTGTCTGACACCGTGGTCGAGCCCTACAATGCCACCCTC
TCCGTCATCAGTTGGTAGAGAATACTGATGAGACCTATTGCATTGACAACGAGGCCCTCTATGATATC
TGCTTCCGCACTCTGAAGCTGACCACACCAACCTACGGGGATCTGAACCACCTTGTCTCAGCCACCATG
AGTGGTGTCAACCCTGCCTCCGTTTCCCTGGCCAGCTCAATGCTGACCTCCGCAAGTTGGCAGTCAAC
ATGGTCCCCTTCCACGTCTCCATTTCTTTATGCCTGGCTTTGCCCTCTCACCAGCCGTGGAAGCCAG
CAGTATCGAGCTCTCACAGTGCCGGAACCTACCCAGCAGGTCTTCGATGCCAAGAACATGATGGCTGCC
TGTGACCCCGCCACGGCCGATACTCACCGTGGCTGTCTTCCGTGGTCGGATGCCATGAAGGAG
GTCGATGAGCAGATGCTTAACGTGCAGAACAAGAACAGCAGCTACTTTGTGGAATGGATCCCCAACAT
GTCAAGACAGCCGCTGTGACATCCCACCTCGTGGCTCAAGATGGCAGTCACCTTCATTGGCAATAGC
ACAGCCATCCAGGAGCTTTCAAGCGCATCTCGGAGCAGTTCACTGCCATGTTCCGCCGGAAGCCCTTC
CTCCACTGGTACACAGGCGAGGGCATGGACGAGATGGAGTTCACCGAGGCTGAGAGCAACATGAACGAC
CTCGTCTCTGAGTATCAGCAGTACCAGGATGCCACCGCAGAAGAGGAGGATTTCGGTGAGGAGGCC
GAAGAGGAGGCCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001293212

Insert Size: 1395 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001293212.1
RefSeq Size:	2772 bp
RefSeq ORF:	1395 bp
Locus ID:	203068
UniProt ID:	P07437
Cytogenetics:	6p21.33
Protein Families:	Druggable Genome
Protein Pathways:	Gap junction, Pathogenic Escherichia coli infection
MW:	52 kDa
Gene Summary:	<p>This gene encodes a beta tubulin protein. This protein forms a dimer with alpha tubulin and acts as a structural component of microtubules. Mutations in this gene cause cortical dysplasia, complex, with other brain malformations 6. Alternative splicing results in multiple splice variants. There are multiple pseudogenes for this gene on chromosomes 1, 6, 7, 8, 9, and 13. [provided by RefSeq, Jun 2014]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a).</p>