

Product datasheet for **SC336283**

CELF3 (NM_001291106) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CELF3 (NM_001291106) Human Untagged Clone
Tag:	Tag Free
Symbol:	CELF3
Synonyms:	BRUNOL1; CAGH4; ERDA4; ETR-1; TNRC4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336283 representing NM_001291106.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGAAGGAGCCGGATGCCATCAAGCTGTTTGTGGGCAGATCCCAGGCATCTGGAGGAGAAGGACCTG
AAGCCCATCTTCGAACAGTTTGGTCGGATCTTTGAGCTGACTGTCATCAAGGACAAGTACACCGGGCTG
CACAAGGGATGTGCTTCTGACATACTGTGCCGGGATTAGCCCTGAAGGCCAGAGCGCCCTGCAC
GAACAGAAGACGCTTCCAGGGATGAACAGGCCGATCCAGGTCAAGCCAGCCGACAGCGAGAGCCGAGGA
GAAGACCGGAAGCTCTTTGTGGGATGCTAGGGAAGCAGCAGACAGATGAGGACGTCGGGAAGATGTTT
GAGCCCTTCGGGACCATCGACGAGTGCCTGTGCTCCGGGGCCAGATGGCACCAGCAAAGGCTGCGCC
TTCGTGAAGTTCCAGACCCAGCTGAGGCCAGGCCGATCAACACCCTTACAGCAGCCGGACCCTG
CCAGGTGCCTCGTCCAGCCTGGTGGTGAAGTTTGTGACTGAGAAGGAGCGAGGTCTCCGCCGATG
CAGCAGGTGGCCACCCAGTTGGGCATGTTAGCCCATCGCCCTCAGTTTGGAGCTACAGCGCCTAC
ACCCAGGCCCTGATGCAGCAGCAGCGGCCCTGGTAGCGGCTCACAGTGCCTACCTAGCCCATGGCC
ACCATGGCTGCCGTGCAGATGCAGCACATGGCTGCCATCAATGCCAATGGCTCATCGCCACCCCATC
ACCCATCCTCAGGAACAGCACCCCTCCTGCCATCGCTGCCACGCCCTGTCTCTGCCATTCGGGCTGCC
CTGGGCGTCAACGGCTACAGCCCGGTGCCACCCAGCCACTGGGCAGCCTGCCCTGATGCTCTGTAT
CCCAACGGGGTTACCCCTACCCAGCCAGAGCCCGCGGCCCGTGGACCCCTGCAGCAGGCCATC
GCGGGGATGCAGCACTACACAGCCTACCCAGCAGCCTACAGCCTGGTGCACCTGCGTTCGCCGAGCCT
CCAGCCCTGGTCGCCAGCAGCCCCACCACCTCAACAGCAGCAGCAGCAGCAGCAGCAGCAACAG
CAGCAGCAGCAAAGAGAAGGCCCTGATGGCTGCAACATCTTATCTACCCTGCCAGGAGTTCACT
GACTCAGAGATCCTCAGATGTTTGTCCCTTTGGCCACGTATCTCAGCCAAAGTCTTTGTTGACCGA
GCCACCAATCAAAGCAAATGTTTGGCTTGTGAGTTTGCACAATCCGGCCAGTGCCAGGCTGCCATC
CAGGCCATGAATGGCTTCCAGATCGGCATGAAGCGCCTCAAAGTCCAGCTAAAGCGGCCTAAGGATGCC
AACCGGCCCTACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001291106

Insert Size: 1395 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001291106.1
RefSeq Size:	5612 bp
RefSeq ORF:	1395 bp
Locus ID:	11189
UniProt ID:	Q5SZQ8
Cytogenetics:	1q21.3
Protein Families:	Transcription Factors
MW:	50.5 kDa
Gene Summary:	<p>Members of the CELF/BRUNOL protein family contain two N-terminal RNA recognition motif (RRM) domains, one C-terminal RRM domain, and a divergent segment of 160-230 aa between the second and third RRM domains. Members of this protein family regulate pre-mRNA alternative splicing and may also be involved in mRNA editing, and translation. Multiple alternatively spliced transcript variants encoding different isoforms have been identified in this gene. [provided by RefSeq, Feb 2010]</p> <p>Transcript Variant: This variant (4) uses an alternate in-frame splice site in the coding region, compared to variant 1. The encoded protein (isoform 4) is shorter than isoform 1.</p>