

Product datasheet for **SC336287**

PRMT5 (NM_001282956) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PRMT5 (NM_001282956) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRMT5
Synonyms:	HRMT1L5; HSL7; IBP72; JBP1; SKB1; SKB1Hs
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336287 representing NM_001282956.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCGGGTCCGAACTCGGGGACGGGAGAGGGCAGACTAGTCATCCCGGAGAAGCAGGGGTTTGATTC
CTCTGCATGCCTGTCTCCATCCGCGTTTCAAGAGGGAGTTTCATTCAGGAACCTGCTAAGAATCGGCC
GGTCCCCAGACAGATCAGACCTACTGCTGT CAGGAAGGACTGGAATACGCTAATTGTGGGAAAGCTT
TCTCCATGGATTCTCCAGACTCAAAGTGGAGAAGATTCGCAGGAACTCCGAGGCGTTGGAGGTGCAG
TTCATCATCACAGGCACCAACCACCTCAGAGAAGGAGTTCTGCTCCTACCTCCAATACCTGGAATAC
TTAAGCCAGAACCCTCCACCTAATGCCTATGAACTCTTTGCCAAGGGCTATGAAGACTATCTGCAG
TCCCCGTT CAGCCACTGATGGACAATCTGGAATCTCAGACATATGAAGTGTGAAAAGGACCCCATC
AAATACTCTCAGTACCAGCAGGCCATCTATAAATGTCTGTAGACCGAGTACCAGAAGAGGAGAAGGAT
ACCAATGTCCAGTACTGATGGTCTGGGAGCAGGACGGGACCCCTGGTGAACGCTTCCCTGCGGGCA
GCCAAGCAGGCCGACCGCGGATAAAGCTGTATGCTGTGGAGAAAACCCAAATGCCGTGGTGACGCTA
GAGAACTGGCAGTTTGAAGAATGGGGAAGCCAAGTGACCGTAGTCTCATCAGACATGAGGGAATGGGTG
GCTCCAGAGAAAGCAGACATCATTGTCAGTGAGCTTCTGGGCTCATTGCTGACAATGAATTGTCGCT
GAGTGCCTGGATGGAGCCAGCACTTCTAAAAGATGATGGTGTGAGCATCCCCGGGGAGTACACTTCC
TTTCTGGCTCCATCTCTTCTCCAAGCTGTACAATGAGGTCCGAGCCTGTAGGGAGAAGGACCGTGAC
CCTGAGGCCAGTTTGAGATGCCTTATGTGGTACGGCTGCACAACCTCCACCAGCTCTCTGCACCCAG
CCCTGTTTACCTTCAGCCATCCAACAGAGATCCTATGATTGACAACAACCGCTATTGCACCTTGAA
TTTCTGTGGAGGTGAACACAGTACTACATGGCTTTGCCGGCTACTTTGAGACTGTGCTTTATCAGGAC
ATCACTCTGAGTATCCGTCAGAGACTCACTCTCTGGGATGTTCTCATGGTTTCCCATCTCTCCCT
ATTAAGCAGCCATAACGGTACGTGAAGGCCAAACCATCTGTGTGCGTTTCTGGCGATGCAGCAATTCC
AAGAAGGTGTGGTATGAGTGGGCTGTGACAGCACCAGTCTGTTCTGCTATTTCATAACCCACAGGCCG
TCATATACCATTGGCCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001282956

Insert Size: 1401 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282956.1
RefSeq Size:	2164 bp
RefSeq ORF:	1401 bp
Locus ID:	10419
UniProt ID:	O14744
Cytogenetics:	14q11.2
Protein Families:	Stem cell - Pluripotency
MW:	53.6 kDa
Gene Summary:	<p>This gene encodes an enzyme that belongs to the methyltransferase family. The encoded protein catalyzes the transfer of methyl groups to the amino acid arginine, in target proteins that include histones, transcriptional elongation factors and the tumor suppressor p53. This gene plays a role in several cellular processes, including transcriptional regulation, and the assembly of small nuclear ribonucleoproteins. A pseudogene of this gene has been defined on chromosome 4. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2015]</p> <p>Transcript Variant: This variant (6) contains multiple differences in the 5' coding region, compared to variant 1, including initiation of translation at an alternate start codon. The encoded isoform (f) is shorter and has a distinct N-terminus, compared to isoform a.</p>