

## Product datasheet for **SC336301**

### ADA2 (NM\_001282227) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ADA2 (NM_001282227) Human Untagged Clone
Tag:	Tag Free
Symbol:	ADA2
Synonyms:	ADGF; CECR1; IDGFL; PAN; SNEDS; VAIHS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC336301 representing NM\_001282227.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGATGCGGGCTGGGGGGCGGCTGGTGTGAACACCAAGGAGGAGCTGGCCAATGAGAGGCTCATGACG
CTCAAAATCGCTGAGATGAAGGAGGCCATGAGGACCCTGATATCCACCAGCATGCACCTTTTCCAG
GCCAAGCATCTCATTGAGAGAAGTCAAGTGTAAATATTCTAAGGATGATGCCAAAAGGGCTGCCTTG
CACCTCCATGACATTGGCATCGTACTATGGACTGGCTGGTGGAGGAATGTACCTACAGGCTCACTGC
CACATCTGTTCACCCCAAGGGGATCATGCAGTTCAGATTTGCTCACCCAACTCCCGTCCATCAGAA
AAATGTTCCAAGTGGATTCTGCTGGAGGATTATCGGAAGCGGGTGCAGAACGTCCTGAGTTTGTGAC
AGCTTGTGAGGAATTTCACTCTGGTGACCCAGCACCCGGAGGTGATTTACACAAACAAAATGTTGTC
TGGTCGAAATTTGAAACCATCTTCTTACCATCTCTGGTCTCATCCATTACGCACCAGTGTTCAGAGAC
TATGTCTTCCGGAGCATGCAGGAGTTCTACGAGGACAACGTGCTCTACATGGAGATCAGAGCCAGGCTG
CTGCCGGTGTATGAGCTCAGTGGAGAGCACCATGACGAAGAGTGGTCAGTGAAGACTTACCAGGAAGTA
GCTCAGAAGTTTGTGGAAACTCACCTGAGTTTATTGGAATCAAATCATTATTTCGGATCACAGATCC
AAAGATGTGGCTGTATCGCAGAATCCATCCGAATGGCCATGGGGCTCCGAATCAAGTTCACCCACGGTG
GTGGCAGGGTTTGACCTGGTGGGCATGAGGACACTGGCCACTCCTTGATGACTACAAGGAAGCTCTG
ATGATCCCGCCAAGGATGGCGTTAAGCTGCCTTACTTCTCCACGCCGGAGAAACAGACTGGCAGGGT
ACTTCCATAGACAGGAACATTCTGGATGCTCTGATGCTGAACACTACCAGAATCGCCATGGATTTGCT
TTGAGCAAACACCCCGCAGTCAGGACTTACTCTGGAAAAGGACATCCCATAGAAGTCTGTCCCATC
TCTAACCAAGTGTGAACTGGTGTCTGACTTGAGGAACCACCTGTAGCCACTGTATGGCCACTGGG
CACCCCATGGTATCAGCTCTGATGACCCAGCTATGTTTGGTGCCAAAGGCTTGTCTATGATTTTCTAT
GAGGTCTTCATGGGCATTGGGGGATGAAGGCTGACCTGAGGACCCTCAAACAGCTGGCCATGAACCTCT
ATCAAGTACAGTACCCTGTTGGAGAGTGAGAAAATACTTTCATGAAATCTGGAAGAAGAGATGGGAT
AAGTTCATAGCAGATGTGGCTACAAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001282227

**Insert Size:** 1410 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001282227.1</a>
<b>RefSeq Size:</b>	4283 bp
<b>RefSeq ORF:</b>	1410 bp
<b>Locus ID:</b>	51816
<b>UniProt ID:</b>	<a href="#">Q9NZK5</a>
<b>Cytogenetics:</b>	22q11.1
<b>Protein Families:</b>	Druggable Genome, Secreted Protein, Transmembrane
<b>MW:</b>	54.4 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of a subfamily of the adenosine deaminase protein family. The encoded protein is one of two adenosine deaminases found in humans, which regulate levels of the signaling molecule, adenosine. The encoded protein is secreted from monocytes undergoing differentiation and may regulate cell proliferation and differentiation. This gene may be responsible for some of the phenotypic features associated with cat eye syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]</p> <p>Transcript Variant: This variant (5) differs in its 5' UTR, and initiates translation at a downstream start codon, compared to variant 3. The encoded isoform (c) has a shorter N-terminus, compared to isoform a. Both variants 5 and 6 encode the same isoform.</p>