

## Product datasheet for **SC336313**

### **YAP1 (NM\_001282099) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	YAP1 (NM_001282099) Human Untagged Clone
Tag:	Tag Free
Symbol:	YAP1
Synonyms:	COB1; YAP; YAP2; YAP65; YKI
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC336313 representing NM\_001282099.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATCCCGGGCAGCAGCCGCCGCTCAACCGGCCCCAGGGCCAAGGGCAGCCGCCTTCGAGCCC
CCGCAGGGGCAGGGCCCGCTCCGGACCCGGGCAACCGGCACCCGCGGCAGCCAGGGCGCCGCGCAG
GCACCCCGCCGGGCATCAGATCGTGCACGTCCGCGGGGACTCGGAGACCGACCTGGAGGCGCTTTC
AACGCCGTCATGAACCCCAAGACGGCCAACGTGCCCCAGACCGTGCCCATGAGGTCCTCGGAAGTGGCC
GACTCCTTCTTCAAGCCGCCGAGCCAAATCCCCTCCGACAGGCCAGTACTGATGCAGGCACTGCA
GGAGCCCTGACTCCACAGCATGTTTCGAGCTCATTCTCTCCAGCTTCTCTGCAGTTGGGAGCTGTTTCT
CCTGGGACTGACCCCACTGGAGTAGTCTCTGGCCAGCAGCTACACCCACAGCTCAGCATCTTCGA
CAGTCTTCTTTGAGATACCTGATGATGTACCTCTGCCAGCAGGTTGGGAGATGGCAAAGACATCTTCT
GGTCAGAGATACTTCTTAAATCACATCGATCAGACAACAACATGGCAGGACCCAGGAAGGCCATGCTG
TCCAGATGAACGTCACAGCCCCACCAGTCCACCAGTGCAGCAGAATATGATGAACTCGGCTTCAGCC
ATGAACCAGAGAATCAGTCAGAGTGCTCCAGTGAACAGGCCACCACCCTGGCTCCCAGAGCCACAG
GGAGGCGTCATGGGTGGCAGCAACTCCAACCAGCAGCAACAGATGCGACTGCAGCAACTGCAGATGGAG
AAGGAGAGGCTGCGGTGAAACAGCAAGAAGTCTTCGGCAGGTGAGGCCACAGGCAATGCGGAATATC
AATCCCAGCACAGCAAATTCCTCAAAATGTCCAGGAGTTAGCCCTGCGTAGCCAGTTACCAACTGGAG
CAGGATGGTGGGACTCAAATCCAGTGTCTTCTCCCGGATGTCTCAGGAATTGAGAACAATGACGACC
AATAGCTCAGATCCTTTCCTTAACAGTGGCACCTATCACTCTCGAGATGAGAGTACAGACAGTGGACTA
AGTATGAGCAGCTACAGTGTCCCTCGAACCCAGATGACTTCTGAACAGTGTGGATGAGATGGATACA
GGTGATACTATCAACAAAGCACCTGCCCTCACAGCAGAACCCTTCCAGACTACCTTGAAGCCATT
CCTGGGACAAATGTGGACCTTGAACACTGGAAGGAGATGGAATGAACATAGAAGGAGAGGAGCTGATG
CCAAGTCTGCAGGAAGCTTTGAGTTCTGACATCCTTAATGACATGGAGTCTGTTTTGGCTGCCACCAAG
CTAGATAAAGAAAGCTTTCTTACATGGTTAAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001282099

**Insert Size:** 1413 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001282099.1](#)

**RefSeq Size:** 5294 bp

**RefSeq ORF:** 1413 bp

**Locus ID:** 10413

**UniProt ID:** [P46937](#)

**Cytogenetics:** 11q22.1

**Protein Families:** Druggable Genome

**MW:** 50.5 kDa

**Gene Summary:** This gene encodes a downstream nuclear effector of the Hippo signaling pathway which is involved in development, growth, repair, and homeostasis. This gene is known to play a role in the development and progression of multiple cancers as a transcriptional regulator of this signaling pathway and may function as a potential target for cancer treatment. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2013]

Transcript Variant: This variant (7) lacks an alternate in-frame exon in the coding region, compared to variant 9. The encoded protein (isoform 7) represents the YAP1-1delta isoform described in Figure 3 of PMID: 22939869. It is shorter, compared to isoform 9. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.