

Product datasheet for **SC336333**

BSDC1 (NM_001300958) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BSDC1 (NM_001300958) Human Untagged Clone
Tag:	Tag Free
Symbol:	BSDC1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336333 representing NM_001300958. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGAAGGGGAGGACGTGGGATGGTGGCGGAGCTGGCTGCAGCAGAGCTACCAAGCAGTCAAAGAG
AAGTCCTCTGAAGCCTTGGAGTTTATGAAGCGGGACCTGACGGAGTTTACCCAGGTGGTGCAGCATGAC
ACGGCCTGTACCATCGCAGCCACGGCCAGCGTGGTCAAGGAGAAGCTGGCTATTGCAGCCTGTTCCCGG
GGCGCTTGCTTCTCTGCCCGTTCTCTATACAGACGGAAGGCTCCTCAGGAGCAACAGAGAAGATGAAG
AAAGGGTTATCTGACTTCTAGGGGTGATCTCAGACACCTTTGCCCTTCGCCAGACAAAACCATCGAC
TGCGATGTATCACCCGTGATGGGCACACCGTCTGGCACAGCTGAGCCCTATGATGGCACCAGGCTCGC
CTCTATAGCCTGCAGTCGGACCCAGCAACCTACTGTAATGAACCAGATGGGCCCCCGGAATTGTTTGAC
GCCTGGCTTTCCAGTTCTGCTTGGAGGAGAAGAAGGGGAGATCTCAGAGCTCCTTGTAGGCAGCCCC
TCCATCCGGGCCCTACACCAAGATGGTTCAGCAGCTGTTTCCATTGAGAAATTCGGCATCGGTAT
TTCTATAAAGTCCATCAGTTAGAGCAGGAGCAGGCCCGGAGGGACGCCCTGAAGCAGCGGGCGGAACAG
AGCATCTCTGAAGAGCCCGGCTGGGAGGAGGAGGAAGAGGAGCTCATGGCATTTCACCCATATCTCCA
AAAGAGGCAAAGGTTCTGTGGCCAAAATTTCTACATTCCTGAAGGAGAACCTGGCCCCAGAGCCCC
TGTGAAGAGAATCTGGTACTTCAGTTGAGCCCCAGCAGAGGTGACTCCATCAGAGAGCAGTGAGAGC
ATCTCCCTCGTGACACAGATCGCCAACCCGGCCACTGCACCTGAGGCACGAGTGCTACCCAAGGACCTG
TCCAAAAGCTGCTAGAGGCATCCTTGGAGGAACAGGGCCTGGCTGTGGATGTGGGTGAGACTGGACCC
TACCCCTATTCACTCCAAGCCCTAACGCCTGCTGGCCACACCGCGGCCAGAGCCCAGGCCTCCA
GCCAGAGTAGAGACTCTGAGGGAGGAGGCGCCACAGACTTACGGGTGTTTGAAGTGAATCGGATAGT
GGGAAGTCTACACCCCTCAACAATGAAAGAAAGGCTCAAGCAGCGGACATCAGTGAGGACTGGGAGAAA
GACTTTGACTTGGACATGACTGAAGAGGAGGTGCAGATGGCACTTTCAAAAGTGGATGCCTCCGGGGAG
GTGAGTGGCCTGGTGGGTGAGAGGGAAGCGAGCCTAATGGTCTGGGTGTGAGAGCTCTCCCCAGCCA
GCCAGCTGTCCCTCAGGAGGGTCCCTGCTCCTGTCTGAGGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001300958
Insert Size:	1425 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001300958.1</u>
RefSeq Size:	5172 bp
RefSeq ORF:	1425 bp
Locus ID:	55108
UniProt ID:	<u>Q9NW68</u>
Cytogenetics:	1p35.1
MW:	51.2 kDa