

## Product datasheet for **SC336337**

### TCP1 theta (CCT8) (NM\_001282909) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	TCP1 theta (CCT8) (NM_001282909) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCT8
Synonyms:	C21orf112; Cctq; D21S246; PRED71
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336337 representing NM\_001282909.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGACCAAAATGGTACAGCATCCTGTGCAAAAATGATTGTAATGGCTTCTCATATGCAAGAGCAAGAA
GTTGGAGATGGCACAACCTTTGTTCTGGTATTTGCTGGAGCTCTCTGGAATTAGCTGAAGAACTTCTG
AGGATTGGCCTGTCAGTTTCAGAGGTCATAGAAGTTATGAAATAGCCTGCAGAAAAGCTCATGAGATT
CTTCTAATTTGGTATGTTGTTCTGCAAAAAACCTTCGAGATATTGATGAAGTCTCATCTCTACTTCGT
ACCTCCATAATGAGTAAACAATATGGAATGAAGTATTCTGGCCAAGCTTATTGCTCAGGCATGCGTA
TCTATTTTTCTGATTCCGGCCATTTCAATGTTGATAACATCAGAGTTTGTAAAATCTGGGCTCTGGT
ATCAGTTCCTCTCAGTATTGCATGGCATGTTTTTAAGAAGGAAACCGAAGGTGATGAACATCTGTC
AAAGATGCAAAAATAGCAGTGTACTCTTGTCTTTTATGATGGCATGATAACAGAACTAAGGGAACAGTG
TTGATAAAGACTGCTGAAGAATTGATGAATTTTAGTAAGGGAGAAGAAAACCTCATGGATGCACAAGTC
AAAGCTATTGCTGATACTGGTGCAAATGTCGTAGTAACAGGTGGCAAAGTGGCAGACATGGCTCTTCAT
TATGCAAATAAATAAATATCATGTTAGTGAGGCTAAACTCAAAATGGGATCTCCGAAGACTTTGTAAA
ACTGTTGGTGCTACAGCTCTTCTAGATTGACACCTCCTGTCCTTGAAGAAATGGGACACTGTGACAGT
GTTTACCTCTCAGAAGTTGGAGATACTCAGGTGGTGGTTTTTAAGCATGAAAAGGAAGTGGCGCCATT
TCTACCATAGTACTTCGAGGCTCTACAGACAATCTGATGGATGACATAGAAAAGGCGAGTAGACGATGGT
GTTAATACTTTCAAAGTCTTACAAGGGATAAACGTCTTGTACCCGGAGGTGGAGCAACAGAAATTGAA
TTAGCCAAACAGATCACATCATATGGAGAGACATGCTCGGACTTGAACAGTATGCTATTAAGAAGTTT
GCTGAGGCATTTGAAGCTATCCCCGCGCACTGGCAGAAAACCTGGAGTTAAGGCAATGAAGTAATC
TCTAAACTTTATGCAGTACATCAAGAAGGAAATAAAAACGTTGGATTAGATATTGAGGCTGAAGCCCT
GCTGTAAGGACATGCTGGAAGCTGGTATTCTAGATACTTACCTGGGAAAATATTGGGCTATCAAACCTC
GCTACTAATGCTGCAGTCACTGTACTTAGAGTGGATCAGATCATCATGGCAAAACAGCTGGTGGGCC
AAGCCTCCAAGTGGGAAGAAAGACTGGGATGATGACCAAATGATGTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN: NM\_001282909

Insert Size: 1428 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001282909.1</a>
<b>RefSeq Size:</b>	1927 bp
<b>RefSeq ORF:</b>	1428 bp
<b>Locus ID:</b>	10694
<b>UniProt ID:</b>	<a href="#">P50990</a>
<b>Cytogenetics:</b>	21q21.3
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	51.6 kDa
<b>Gene Summary:</b>	<p>This gene encodes the theta subunit of the CCT chaperonin, which is abundant in the eukaryotic cytosol and may be involved in the transport and assembly of newly synthesized proteins. Alternative splicing results in multiple transcript variants of this gene. A pseudogene related to this gene is located on chromosome 1. [provided by RefSeq, Sep 2013]</p> <p>Transcript Variant: This variant (4) lacks an alternate exon in the 5' coding region and initiates translation at an alternate start codon, compared to variant 1. It encodes isoform 4, which has a shorter and distinct N-terminus, compared to isoform 1.</p>