

Product datasheet for **SC336339**

TMEM161B (NM_001289008) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM161B (NM_001289008) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM161B
Synonyms:	FLB3342; PRO1313
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336339 representing NM_001289008.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGTGATGGCCAGTGTATGCAGAAGATTATACCTCACTATTCTTGTCTCGATGGCTACTCTGTAAT
GGCAGTTTGAGGTGGTATCAACATCCTACAGAAGAATTAAAGATTCTTGCAGGGAACAACAAAAA
GGGAAAAACAAAAAGATAGGAAATATAATGGTCACATTGAAAGTAAGCCATTAACCATTCCAAAGGAT
ATTGACCTTCATCTAGAAACAAAGTCAGTTACAGAAGTGATACTTTAGCATTGCATTACTTTCCAGAA
TACCAGTGGCTGGTGGATTTACAGTGGCTGCTACAGTTGTGTATCTAGTAACTGAAGTCTACTACAAT
TTTATGAAGCCTACACAGGAAATGAATATCAGCTTAGTCTGGTGCCTACTTGTCTTTGCAATC
AAAGTTCTATTTTCATTAACACACTATTTTAAAGTAGAAGATGGTGGTGAAAGATCTGTTTGTGTC
ACCTTTGGATTTTTTTCTTTGTCAAAGCAATGGCAGTGTGATTGTAACAGAAAATTATCTGGAATTT
GGACTTGAACAGGGTTTACAAATTTTTCAGACAGTGCAGTGCAGTTTCTGAAAAGCAAGGTTAGAA
TCTCAGAGTCTGTTTCAAACTTACTTTCAAATTTTCTGGCTATTTTCTGTTCACTTCTGGGGCT
TTTTTGACATTTCTGGATTACGACTGGCTCAAATGCATCTGGATGCCCTGAATTTGGCAACAGAAAAA
ATTACACAACTTTACTTCATATCAACTTCTTGGCACCTTTATTTATGGTCTGCTCTGGGTAAACCA
ATCACCAAGACTACATTATGAACCCACCACTGGGCAAGAAAGTATCCCTTTAATGACAGAAGCCACA
TTCGATACTCTGCGACTCTGGTTAATAATCCTGCTGTGTGCTTTGCGGTTGGCCATGATGCGTAGTCAC
CTGCAAGCTTATTTAAATTTAGCCCAAAATGTGTGGATCAGATGAAGAAAGAAGCGGGGCGAATAAGC
ACGGTTGAGCTACAGAAAATGGTGGCTCGAGTCTTTTATTATCTTTGTGTCATTGCACTGCAGTATGTG
GCGCCTCTGGTAATGCTGCTTACACAACTCTGCTTTTGAACAACTAGGTAATCATTCTGGGGTATT
TATCCAGAATCTATCTACCTTACCAGTGGATAATAGTCTACTGTCCAATTCTGTTACTCTGAATTA
CCATCAGCTGAAGGGAATGAAGGTAAGTGTACACAAATAACAGTGGCACTGAGCAGCTTAAAAAT
ATTTTACTCCTCTTCTTTTTCGAGGACTTCTGTCTTTCTGACCTGGTGGATTGCTGCTTGCCTCTTT
TCTACAAGCCTTTTGGGCTTTTCTATCACCAGTATCTGACTGTGGCATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001289008

Insert Size: 1431 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001289008.1

RefSeq Size:	2766 bp
RefSeq ORF:	1431 bp
Locus ID:	153396
Cytogenetics:	5q14.3
Protein Families:	Transmembrane
MW:	54.4 kDa