

Product datasheet for **SC336355**

PCYOX1L (NM_001301055) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PCYOX1L (NM_001301055) Human Untagged Clone
Tag:	Tag Free
Symbol:	PCYOX1L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336355 representing NM_001301055. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCCGCCGGGCAAAATCGGTGCGGGGAAGGACGCGCGGTGGTTGGGGCTGGGATTGGGGCTCTGCT
GTGGCCATTTTCTCCAGCAGCACTTTGGACCTCGGGTGCAGATCGACGTGTACGAGAAGGGAACCGTG
GGTGGCCGCTTGGCCACCATCTCAGTCAACAAGCAGCACTATGAGAGCGGGGCTGCCTCCTTCCACTCC
CTGAGCCTGCACATGCAGGACTTCGTCAAGCTGCTGGGCTGAGGCACCGCGCGAGGTGGTGGGCAGG
AGCGCCATCTTCGGCGGGGAGCACTTCATGCTGGAGGAGACTGACTGGTACCTGCTGAACCTTTCCGC
CTCTGGTGGCACTATGGCATCAGCTTCTGAGGCTGCAGATGTGGGTGGAGGAGGTCATGGAGAAGTTC
ATGAGGATCTATAAGTACCAGGCCACGGCTATGCCTTCTCGGGTGTGGAGGAGCTGCTCTACTCACTG
GGGAGTCCACCTTTGTTAACATGACCCAGCACTCTGTGGCTGAGTCCCTGCTGCAGGTGGGCGTCACG
CAGCGCTTATTGATGATGTCGTTTCTGCTGTCCTGCGGGCCAGCTATGGCCAGTCAGCAGCGATGCC
GCCTTTGAGGAGCCATGTCACTAGCCGGGCCCCAAGGCAGCCTGTGGTCTGTGGAAGGAGGCAATAAG
CTGTTTGTTCGGTTTCTGAAGCTCACCAAGGCCAATGTGATCCATGCCACAGTGACCTCTGTGACC
CTGCACAGCACAGAGGGAAAGCCCTGTACCAGGTGGCGTATGAGAATGAGGTAGGCAACAGCTCTGAC
TTCTATGACATCGTGGTTCATCGCCACCCCTGCACCTGGACAACAGCAGCAGCAACTTAACCTTTGCA
GGCTTCCACCCGCCATTGATGACGTGCAGGGCTTTTCCAGCCACCGTCTCTCTTGGTCCACGGC
TACCTCAACTCGTCTACTTCGGTTTCCAGACCCCTAAGCTTTTCCCTTTGCCAACATCCTTACCACA
GATTTCCCGAGCTTCTTCTGCACTCTGGACAACATCTGCCCTGTCAACATCTCTGCCAGTTCGGCGGA
AAGCAGCCCCAGGAGGAGCTGTTTGGCGAGTCCAGTCCCCAAGCCCTCTTTCGGACCCAGCTAAAG
ACCCTGTTCCGTTCTATTACTCAGTGCAGACAGCTGAGTGGCAGGCCATCCCCCTATGGCTCCCGC
CCCAGCTCCCGAGGTTTGCACCTCATGACCAGCTTTTACCTCAATGCCCTGGAGTGGGCGGCCAGC
TCCGTGGAGGTGATGGCCGTGGCTGCCAAGAATGTGGCCTTGCTGGCTTACAACCGCTGGTACCAGGAC
CTAGACAAGATTGATCAAAAAGATTTGATGCACAAGGTCAAGACTGAAGTGATGA
ACGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001301055

Insert Size: 1434 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001301055.1](#)

RefSeq Size: 2583 bp

RefSeq ORF: 1434 bp

Locus ID: 78991

Cytogenetics: 5q32

MW: 53.4 kDa

Gene Summary:

Probable oxidoreductase.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) represents the longest transcript. It has an alternate splice site in the 5' region, which results in an alternate translation initiation site, compared to variant 1. The resulting isoform (2) is shorter and has a distinct N-terminus, compared to isoform 1.