

Product datasheet for **SC336359**

IFIT1 (NM_001270927) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IFIT1 (NM_001270927) Human Untagged Clone
Tag:	Tag Free
Symbol:	IFIT1
Synonyms:	C56; G10P1; IFI-56; IFI-56K; IFI56; IFIT-1; IFNAI1; ISG56; P56; RNM561
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336359 representing NM_001270927.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGTACAAATGGTGATGATCATCAGGTCAAGGATAGTCTGGAGCAATTGAGATGTCACTTTACATGG
GAGTTATCCATTGATGACGATGAAATGCCTGATTTAGAAAACAGAGTCTTGGATCAGATTGAATTCCTA
GACACCAAATACAGTGTGGGAATACACAACCTACTAGCCTATGTGAAACACCTGAAAGGCCAGAATGAG
GAAGCCCTGAAGAGCTTAAAAGAAGCTGAAAACCTAATGCAGGAAGAACATGACAACCAAGCAAATGTG
AGGAGTCTGGTACCTGGGGCAACTTTGCCTGGATGTATTACCACATGGGCAGACTGGCAGAAGCCAG
ACTTACCTGGACAAGGTGGAGAACATTTGCAAGAAGCTTTCAAATCCCTCCGCTATAGAATGGAGTGT
CCAGAAATAGACTGTGAGGAAGGATGGGCCTTGTGAAGTGTGGAGGAAAAAATTATGAACGGGCCAAG
GCCTGCTTTGAAAAGTGCTTGAAGTGGACCCTGAAAACCCTGAATCCAGCGCTGGGTATGCGATCTCT
GCCTATCGCCTGGATGGCTTTAAATTAGCCACAAAAAATCACAAGCCATTTTCTTTGCTTCCCCTAAGG
CAGGCTGTCCGCTTAAATCCAGACAATGGATATATTAAGGTTCTCCTTGCCTGAAGCTTCAGGATGAA
GGACAGGAAGCTGAAGGAGAAAAGTACATTGAAGAAGCTCTAGCCAACATGTCCTCACAGACCTATGTC
TTTCGATATGCAGCCAAGTTTTACCGAAGAAAAGGCTCTGTGGATAAAGCTCTTGAGTTATAAAAAAG
GCCTTGCAGGAAACCCCACTTCTGTCTTACTGCATCACCAGATAGGGCTTTGCTACAAGGCACAAATG
ATCCAAATCAAGGAGGCTACAAAAGGCCAGCCTAGAGGGCAGAACAGAGAAAAGCTAGACAAAATGATA
AGATCAGCCATATTTTCAATTTGAATCTGCAAGTGGAAAAAAGCCACATTTGAGGTGGCTCATCTAGAC
CTGGCAAGAAATGTATATAGAAGCAGGCAATCACAGAAAAGCTGAAGAGAAATTTCAAATTTGTTATGC
ATGAAACCAGTGGTAGAAGAAAACAATGCAAGACATACATTTCCACTATGGTCGGTTTCAGGAATTTCAA
AAGAAATCTGACGTCAATGCAATTAATCCATTATTTAAAAGCTATAAAAAATAGAACAGGCATCATTAAAC
AGGGATAAAAAGTATCAATTCCTTTGAAGAAATTGGTTTTAAGGAAACTTCGGAGAAAAGGCATTAGATCTG
GAAAGCTTGAGCCTCCTTGGGTTCTGCTACAAATTGGAAGGAAATATGAATGAAGCCCTGGAGTACTAT
GAGCGGGCCCTGAGACTGGCTGCTGACTTTGAGAACTCTGTGAGACAAGGCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN: NM_001270927

Insert Size: 1437 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270927.1</u>
RefSeq Size:	4493 bp
RefSeq ORF:	1437 bp
Locus ID:	3434
Cytogenetics:	10q23.31
MW:	55.4 kDa
Gene Summary:	<p>This gene encodes a protein containing tetratricopeptide repeats that was originally identified as induced upon treatment with interferon. The encoded protein may inhibit viral replication and translational initiation. This gene is located in a cluster on chromosome 10 with five other closely related genes. There is a pseudogene for this gene on chromosome 13. Alternatively spliced transcript variants encoding multiple isoforms have been observed. [provided by RefSeq, Aug 2012]</p> <p>Transcript Variant: This variant (2) contains an alternate exon in the 5' coding region, and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (2) has a distinct N-terminus, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>