

Product datasheet for **SC336366**

C9orf43 (NM_001278630) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C9orf43 (NM_001278630) Human Untagged Clone
Tag:	Tag Free
Symbol:	C9orf43
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC336366 representing NM_001278630. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAAGCACTAGAATGCTGCAGGGTTGGCCTTTGGCCTAAACCATTTCTAGCTATGGACTTGCCAGAT
GAGAGCCAGTGGGACGAAACCACCTGTGGCTTGGCTGTTGTGCAGCACCCACAATGCTGGCAACTATC
CGCCGATTGAGAGGGGCCATCCTCGAATCCTCGGCTCATCTGCAAAATCCCTGGATGCTGAAGAT
AACTCCCAGTGCTCACCGTGGTAGACATCTTAGATCCGGCTTTCAGCTCATCATTTACCAGAATGT
ACCTTTACTAAGGCCATTCTTTATTGTCTCAGAGTTCAAAGTTTACTCCAAATTTCCAGTGTGAAT
CCGAAGGTTTACCTGACAAAAGTTTGATCAACTGTACTAACAGACTTCCCAAATTTCCAGTGTGAAT
TTGAATGAGACGCAACTTCCCTGCCCTGAAGATGTTAGAAATATGGTTGATTGTGGATCCAGAAGAA
ACAGAGATACATGTGAGCCAGCATGGGAAGAAGAAAAGAAAGAACTCGGCAGTGAAGCAAGTCATTT
CTGGGTCTCTCTGGAAATCAGTCCGAGGAACACGAGTAGGAACACCAGGGATGATCGTGCCTCCCCCA
ACCCAGTGCAATTGTCTGAACAATTCAGTTCAGATTTCTACCTCTCTGGGCTCAATCCGAAGCGTTA
CCTCAGGATCTACTGAAGGAACTTTGGCAGGTGGAAAGCAAACCATGCTCTGTCCAGAGATGAAGATA
AAATTGGCCATGATGAAAAGAATCTTCCCTTGGAAAAGAACCGACCTGACAGTGTGATTTCTTCTAAG
ATGTTTCTATCTATACCCGCTCACCTGGAAAGACCAGCACTGCGATATCCTGAACGTTTGAAGAAA
TTACATAACCTGAAGACAGAAGTTACAGGAAACAGCAGCAGCGGAGCAGCAGCAGCAGCAGCAACAG
AAGAAGGTGAAAACACCTATTAAGAAACAGGAGGCTAAAAAGAAAGCCAAGAGTGATCCAGGGATCCAG
AGCACTTACATAAACATCCAGTTACCACCGTTTCATGACCGTCTCTATGGTTACAGAACTCTGCCAGGT
CAGAACAGTGACATGAAGCAGCAGCAGCAGATGGAAAAGGAACCACTTCGAAACAGGATTCACCGGAG
AGACCAAGATGAACTACTATGACCATGCGGATTTCCACCACAGTGTAAGAAAGTCTGAATTGTATGAA
ACAGAACCCTAACAAAGGACATTAGTGCTCCAGTGGACGCTGTGCCAGAAGCCAGGCTGCCAGGCAA
AAGAAGATCTCCTTTAACTTTTCAGAAATATGGCTAGCACAGGCTGGAACCTCTGAGCTCAAACACTT
AGGATTTCTCAGGACACTGATGATGAGGATGAGGAGGACCAGTCTCTGGGCAGAGTGA
ACGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001278630

Insert Size: 1440 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278630.1](#)

RefSeq Size: 1948 bp

RefSeq ORF: 1440 bp

Locus ID: 257169

Cytogenetics: 9q32

MW: 54.3 kDa