

Product datasheet for **SC336370**

HARS1 (NM_001289094) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HARS1 (NM_001289094) Human Untagged Clone
Tag:	Tag Free
Symbol:	HARS1
Synonyms:	CMT2W; HARS; HRS; USH3B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336370 representing NM_001289094.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATCGAGGAGGAGGTGGCGAAACTCCTGAAACTGAAGGCACAGCTGGGTCTGATGAAAGCAAACAG
AAATTTGTGCTCAAACCCCAAGGGCACAAGAGACTATAGTCCCAGGAGATGGCAGTTGCGGAGAAG
GTGTTTGACGTAATCATCCGTTGCTTCAAGCGCCACGGTGCAGAAGTCATTGATACACCTGTATTTGAA
CTAAAGGAAACACTGATGGGAAAGTATGGGAAAGACTCCAAGCTTATCTATGACCTGAAGGACCAGGGC
GGGAGACTCCTGTCCCTTCGCTATGACCTCACTGTTTCTTTGCTCGGTATTTGGCAATGAATAAATG
ACCAACATTAACGCTACCACATAGCAAAGGTATATCGGCGGGATAACCCAGCCATGACCCGTGGCCGA
TACCGGAATTCTACCAGTGTATTTGACATTGCTGGGAACTTTGATCCCATGATCCCTGATGCAGAG
TGCTGAAGATCATGTGCGAGATCCTGAGTTCACCTCAGATAGGGGACTTCTGGTCAAGGTAACGAT
CGACGCATTCTAGATGGGATGTTTGTATCTGTGGTGTCTGACAGCAAGTCCGTACCATCTGCTCC
TCAGTAGACAAGCTGGACAAGGTGCTCGGGAAGAGGTGAAGAATGAGATGGTGGGAGAGAAGGGCCTT
GCACCTGAGGTGGCTGACCGCATTGGGGACTATGTCCAGCAACATGGTGGGGTATCCCTGGTGGAAACAG
CTGCTCCAGGATCCTAAACTATCCAAAACAAGCAGGCCTTGGAGGGCCTGGGAGACCTGAAGTTGCTC
TTTGAGTACCTGACCCTATTTGGCATTGATGACAAAATCTCCTTTGACCTGAGCCTTGCTCGAGGGCTG
GATTACTACACTGGGGTATCTATGAGGCAGTGTCTACAGACCCAGCCAGGCAGGGGAAGAGCCC
CTGGGTGTGGGCAGTGTGGCTGCTGGAGGACGCTATGATGGGCTAGTGGGCATGTTGACCCCAAAGGG
CGCAAGGTGCCATGTGTGGGCTCAGCATTGGGGTGGAGCGGATTTTCTCCATCGTGGAACAGAGACTA
GAGGCTTTGGAGGAGAAGATACGGACCAGGAGACACAGGTGCTTGTGGCATCTGCACAGAAGAAGCTG
CTAGAGGAAAGACTAAAGCTTGTCTCAGAAGTGGGATGCTGGGATCAAGGCTGAGCTGAGTGTACAAG
AAGAACCCAAAGCTACTGAACCAGTTACAGTACTGTGAGGAGGCAGGCATCCCACTGGTGGCTATCATC
GGCGAGCAGGAACTCAAGGATGGGGTATCAAGCTCCGTTCAAGTACGAGCAGGGAAGAGGTGGATGTC
CGAAGAGAAGACCTTGTGGAGGAAATCAAAGGAGAACAGGCCAGCCCCTCTGCATCTGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN: NM_001289094

Insert Size: 1443 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001289094.1
RefSeq Size:	2219 bp
RefSeq ORF:	1443 bp
Locus ID:	3035
Cytogenetics:	5q31.3
Protein Pathways:	Aminoacyl-tRNA biosynthesis
MW:	54.2 kDa
Gene Summary:	<p>Aminoacyl-tRNA synthetases are a class of enzymes that charge tRNAs with their cognate amino acids. The protein encoded by this gene is a cytoplasmic enzyme which belongs to the class II family of aminoacyl-tRNA synthetases. The enzyme is responsible for the synthesis of histidyl-transfer RNA, which is essential for the incorporation of histidine into proteins. The gene is located in a head-to-head orientation with HARSL on chromosome five, where the homologous genes share a bidirectional promoter. The gene product is a frequent target of autoantibodies in the human autoimmune disease polymyositis/dermatomyositis. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012]</p> <p>Transcript Variant: This variant (7) uses an alternate in-frame splice site in the 5' coding region, compared to variant 1. The encoded protein (isoform 7) lacks an internal region and is shorter, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>