

Product datasheet for **SC336381**

LILRB2 (NM_001278404) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LILRB2 (NM_001278404) Human Untagged Clone
Tag:	Tag Free
Symbol:	LILRB2
Synonyms:	CD85D; ILT-4; ILT4; LIR-2; LIR2; MIR-10; MIR10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336381 representing NM_001278404.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACAGGAGCCTACCCAAAACCCACCCTCTCAGCCCAGCCCAGCCCTGTGGTGACCTCAGGAGGAAGG
GTGACCCCTCCAGTGTGAGTCACAGGTGGCATTGGCGGCTTCATTCTGTGTAAGGAAGGAGAAGATGAA
CACCCACAATGCCTGAACTCCCAGCCCCATGCCGTGGGTGTCGCCGCGCCATCTTCTCCGTGGGCCCC
GTGAGCCCGAATCGCAGGTGGTCGCACAGGTGCTATGGTTATGACTTGAAGTCTCCCTATGTGTGGTCT
TCACCCAGTGATCTCCTGGAGCTCCTGGTCCCAGGTGTTTCTAAGAAGCCATCACTCTCAGTGCAGCCG
GGTCTGTGATGGCCCTGGGAAAGCCTGACCCTCCAGTGTGTCTGTGATGTCGGCTATGACAGATTT
GTTCTGTACAAGGAGGGGAACGTGACCTTCGCCAGCTCCCTGGCCGGCAGCCCCAGGCTGGGCTCTCC
CAGGCCAACTTACCCTGGGCCCTGTGAGCCGCTCCTACGGGGCCAGTACAGATGCTACGGTGCACAC
AACCTCTCTGAGTGTGCTGGCCCCCAGCGACCCCTGGACATCCTGATCACAGGACAGATCCGTGGC
ACACCTTCATCTCAGTGCAGCCAGGCCCCACAGTGGCCTCAGGAGAGAAGCTGACCCTGCTGTGTCAG
TCATGGCCGCGAGTCCACACTTTCTTCTGACCAAGGCGGGAGCAGCTGATGCCCCACTCCGTCTAAGA
TCAATACACGAATATCCTAAGTACCAGGCTGAATTCCCCATGAGTCTGTGACCTCAGCCCACGCGGGG
ACCTACAGGTGCTACGGCTCACTCAACTCCGACCCCTACCTGCTGTCTCACCCCAGTGAGCCCTGGAG
CTCGTGGTCTCAGGACCCTCCATGGGTTCCAGCCCCCACCACCCGGTCCCATCTCCACACCTGCAGGC
CCTGAGGACCAGCCCCTACCCCCACTGGGTCGGATCCCCAAAGTGGTCTGGGAAGGCACCTGGGGGT
GTGATCGGCATCTGGTGGCCGTCGCTACTGCTCCTCCTCCTCCTCCTCCTCCTCCTCCTCCTCCTCCGA
CATCGAGTCCAGGCCAACACTGGACATCGACCCAGAGAAAGGCTGATTTCCAACATCCTGCAGGGGCT
GTGGGCCAGAGCCCACAGACAGAGGCCCTGCAGTGGAGTCCAGCCAGCTGCCGACGCCAGGAAGAA
AACCTCTATGCTGCCGTGAAGGACACACAGCCTGAAGATGGGGTGGAGATGGACACTCGGCTGCTGCA
TCTGAAGCCCCCAGGATGTGACCTACGCCAGCTGCACAGTTCAGCCCTCAGACGGAAGGCAACTGAG
CCTCCTCCATCCCAGAAAGGGAACCTCCAGTGCAGCCAGCATCTACGCCACCCTGGCCATCCACTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001278404

Insert Size:	1449 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001278404.2</u>
RefSeq Size:	2789 bp
RefSeq ORF:	1449 bp
Locus ID:	10288
Cytogenetics:	19q13.42
Protein Families:	Transmembrane
MW:	52 kDa
Gene Summary:	<p>This gene is a member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at chromosomal region 19q13.4. The encoded protein belongs to the subfamily B class of LIR receptors which contain two or four extracellular immunoglobulin domains, a transmembrane domain, and two to four cytoplasmic immunoreceptor tyrosine-based inhibitory motifs (ITIMs). The receptor is expressed on immune cells where it binds to MHC class I molecules on antigen-presenting cells and transduces a negative signal that inhibits stimulation of an immune response. It is thought to control inflammatory responses and cytotoxicity to help focus the immune response and limit autoreactivity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (4) lacks a portion of the 5' coding region, and uses a downstream in-frame start codon, compared to variant 1. The encoded isoform (3) has a shorter N-terminus, compared to isoform 1.</p>