

Product datasheet for **SC336382**

WDR86 (NM_001284260) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WDR86 (NM_001284260) Human Untagged Clone
Tag:	Tag Free
Symbol:	WDR86
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336382 representing NM_001284260.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGGGGCGGGTTCGGCCCTGAGGGTCTGCGCCGACCACCGGGGGCATCAACTGGCTGAGCCTG
AGCCCGACGGGACGGCCTGCTGACGGGCAGCGAGGACGGCACGGCCCGGCTCTGGAGCACCGCGGAC
GGCCAGTGCTGCGCGCTCCTGCAAGGACATGAAAGCTATGTGACCTTCTGCCAGCTGGAGGATGAGGCT
GCCTTCACATGCAGCGCCGACTGCACCATCAGGAGGTGGGACGTGCTGACCGGGCAGTGTCTGCAGGTG
TACCGAGGACACAGTCCATCGTGAACAGGATCCTGGTTGCCAACACCAGCTCTTCAGCAGCTCCTAT
GACCGGACAGCTCGGGTCTGGAGTGTGGACAAGGGGCAGATGTCCGGGAGTTCGGGGCCACCGCAAC
TGCGTGTGACCTAGCCTACTCTGCCCGTGGGACCTCCCAGCACTCCTGCGCGGAGGAGGCCCGG
GCCGGGGGCTTCTGGTGACCGGCAGCACAGATGGCACAGCCAAGGTGTGGCAGGTGGCCAGCGGCTGC
TGCCACCAGACGCTGCGGGCCACACGGGTGCAGTGTGCTGTGCTAGTGTAGACACGCCCGGCCACAG
GCCTTCACAGGCAGCACCGGACGCCACCATCCGTGCTGGGACATCCTGAGTGGGAGCAGCTGCGGGTG
TTCCGGGAGCACCGGGGCTCCGTCTGTCTGGAGAGGGAGAGGTCTGGGGGGCTTCAGAGCTGCCCC
AGCTGCCGACCTTCCCTGCTCCTCCTCGCAGCTGGTGAACCGACTCGTGTACTCTGGCAGCGCGGACAG
GACCGTCAAGTGTCTGGTGGCAGACACAGGGGAGTGTGTGCGCACGTTACGCGCCACAGACGCAACGT
GAGCGCCCTCAAGTACCACGCGGGCACCTTGTTCACGGGCAGCGGGGACGCTTGCGCCGGGGCTTCGA
CGCGCAGTCTGGAGAGCTGCGGAGGGTGTTCGGGGCCACACATTCATCATCAACTGCATCCAGGTGCA
CGGCCAGTGTCTACACCGCTCGCACGACGGCGCCCTGCGCCTCTGGGACGTGCGCGGGCTCCGAGG
TGCCCCGCGCCCTCCGCCATGCGCAGCCTCTCGCGGCTCTCAGCAACAAGGTGGGCTGCGCCCG
CGCGCCCTGCAGCGGCTGATCCCGCGGGGCCCTGCAGACGCCAGCCAGACACCCAGCGGCTCCC
AGAGCGCCCGCCCTGCTACCCGCGGTGGTGGCGCCCGATGGCCGGCAGGGGGCAGGAGCGAGGAAGC
CCGGGCGGGAGGAGGCCGTCGAGGCGTCTGGTTTTTCTTTGGTGGCCAGGAGCGCTGGGAGCGGG
AGTGCTCGCCCTGGGACCGCCCTTTTCCCTTTAGGGTGGCTCCTGTCTCCCTCCCATCCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001284260

Insert Size: 1449 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001284260.1
RefSeq Size:	2120 bp
RefSeq ORF:	1449 bp
Locus ID:	349136
UniProt ID:	Q86T14
Cytogenetics:	7q36.1
MW:	52.2 kDa