

Product datasheet for **SC336395**

Noelin (OLFM1) (NM_001282611) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Noelin (OLFM1) (NM_001282611) Human Untagged Clone
Tag:	Tag Free
Symbol:	OLFM1
Synonyms:	AMY; NOE1; NOELIN1; OlfA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336395 representing NM_001282611.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCGGTGCCGCTGCTCAAGATCGGGGTCGTGCTGAGCACCATGGCCATGATCACTAACTGGATGTCC
CAGACGCTGCCCTCGCTGGTGGGCTCAACACCACCAAGCTCTCGGCGGCCGGCGGGGACGCTGGAC
CGCAGCACCGGCGTGTGCCACCAACCTGAGGAGAGCTGGCAGGTGTACAGCTCTGCCAGGACAGC
GAGGGCAGGTGTATCTGCACAGTGGTCCACAGCAGACCATGTGTTACGGGATGCCCGCACAAAA
CAGCTGAGGCAGCTACTGGAGAAGTGCAGAACATGTCTCAATCCATAGAGGTCTTGACAGGCGGACC
CAGAGAGACTTGCAGTACGTGGAGAAGATGGAGAACAAATGAAAGGACTGGAGTCCAAGTTCAAACAG
GTGGAGGAGAGTCATAAGCAACACCTGGCCAGGCAGTTTAAAGGCGATAAAAAGCGAAAATGGATGACTT
AGGCCTTTGATACCTGTGTTGGAAGAGTACAAGGCCGATGCCAAATTGGTATTGCAGTTTAAAGAGGAG
GTCCAGAACTGACGTCAGTGTAAACGAGCTGCAAGAGGAAATGGCGCCTATGACTACGATGACTT
CAGAGCAGAGTGTCCAATCTTGAAGAAAGCTCCGTGCATGCATGCAAAAAGTCTTGGCGGAAGTTG
ACGGGCATCAGTGACCCCGTACTGTCAAGACCTCCGGCTCGAGGTTCCGATCCTGGATGACAGACCCT
CTCGCCCTGAAGGCGATAACCGGGTGTGTACATGGACGGCTATCACAACAACCGCTTCGTACGTGAG
TACAAGTCCATGGTTGACTTCAAGAACCGGACAATTTCACTCCACCGTCTCCCCACCCCTGGTTCG
GGCAGCGGGCAGGTGGTCTACAACGGTTCTATCTACTTCAACAAGTTCAGAGCCACATCATCATCAGG
TTTGACCTGAAGACAGAGACCATCTCAAGACCCGAGCTGGACTATGCCGGTTACAACAACATGTAC
CACTACGCTGGGGTGGCCACTCGGACATCGACCTCATGGTGGACGAGAGCGGGCTGTGGCCGTGTAC
GCCACCAACCAGAACGCTGGCAACATCGTGGTCACTAGGCTGGACCCCGTGTCCCTGCAGACCCTGCAG
ACCTGGAACACGAGCTACCCCAAGCGCAGCGCCGGGAGGCCTTCATCATCTGCGGCACGCTGTACGTC
ACCAACGCTACTCAGGGGTACCAAGTCCACTATGCATACCAGACCAATGCCTCCACCTATGAATAC
ATCGACATCCCATTCCAGAACAATACTCCACATCTCCATGCTGGACTACAACCCCAAGGACCGGGCC
CTGTATGCCTGGAACAACGGCCACCAGATCCTCTACAACGTGACCCTCTCCACGTATCCGCTCCGAC
GAGTTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001282611

Insert Size:	1458 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282611.1
RefSeq Size:	2791 bp
RefSeq ORF:	1458 bp
Locus ID:	10439
UniProt ID:	Q99784
Cytogenetics:	9q34.3
Protein Families:	Druggable Genome
MW:	55.3 kDa
Gene Summary:	<p>This gene product shares extensive sequence similarity with the rat neuronal olfactomedin-related ER localized protein. While the exact function of the encoded protein is not known, its abundant expression in brain suggests that it may have an essential role in nerve tissue. Several alternatively spliced transcripts encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR and 5' coding region which results in the use of an alternate start codon, compared to variant 1. The encoded isoform (4) has distinct N-terminus, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>