

Product datasheet for **SC336398**

RBM46 (NM_001277171) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RBM46 (NM_001277171) Human Untagged Clone
Tag:	Tag Free
Symbol:	RBM46
Synonyms:	CT68
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336398 representing NM_001277171.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAATGAAGAAAATATAGATGGAACAAATGGATGCAGTAAAGTTCGAACTGGTATTGCAATGAAGCA
GCATTACTTGGTTGATGGAAAAGACTGGTTACAACATGGTTCAGGAAAATGGACAAAGGAAATTTGGC
GGTCCTCCTCCAGGTTGGGAAGGTCACCTCCACCTAGAGGCTGTGAAGTTTTGTAGGAAAAATACCT
CGTGATATGTATGAAGATGAGTTAGTTCCTGTATTTGAAAGAGCTGGGAAGATATATGAATTTGCACTT
ATGATGGAATTTAGTGGTAAAAATCGAGGTTATGCTTTTGTGATGTACACTACAAAAGAAGAAGCCCAA
TTAGCCATCAGAATTCCTAATAATTATGAAATTCGACCAGGGAAGTTTATTGGTGTGTGTGAAGCCTG
GATAATTGTAGATTATTTATTGGAGCTATCCCAAGGAAAAGAAGAAAGAAGAAATTTAGATGAAATG
AAGAAAGTTACAGAAGGAGTTGTAGATGCATTGTTTATCCAAGTCAACTGATAAGACCAAAAAATCGT
GGTTTTGCATTTGTGGAATATGAATCTCACAGAGCTGCTGCTATGGCAAGGAGGAAACTAATCCAGGA
ACATCCAACATATGGGGCCACACCATTCAGGTAGATTGGGCTGACCCAGAGAAAGAGGTGGATGAGGAA
ACCATGCAGAGAGTTAAAGTTCTTTATGTAAGAAATTTAATGATCTCAACTACAGAGGAAACAATTTAA
GCAGAATTCATAAATTTAAGCCTGGTGCAGTTGAACGGGTAAAGAACTTAGAGATTATGCTTTTGT
CACTTTTTCAACCGAGAAGATGCAGTGGCTGCCATGTCTGTTATGAATGGAAAATGCATTGATGGAGCA
AGTATTGAGGTAACACTAGCTAAACCAGTAAATAAAGAAAACACTGGAGACAGCATCTTAATGGTCAG
ATTAGTCCAAATCTGAAAATCTGATTGTGTTTGTCAACAAGAAGAGAGCCACCCAAAAACTTAGGC
AAGCTGCCAACTCTTCTGCTCGTCTCAATGGTCAGCATAGCCCAAGTCCGCCTGAAGTTGAAAGATGC
ACTTACCCTTTTTATCCTGGAACAAAGCTTACTCCAATTAGTATGTATTCTTTAAATCCAATCATTTT
AATTCTGCAGTAATGCATTTGGATTATTACTGCAACAAAAATAACTGGGCACCACCAAGATATTATTTA
TATTCAACAACAAGTCAAGATGGGAAAGTACTCTTGGTGTATAAGATAGTTATTCCTGCTATTGCAAT
GGATCCCAGAGTTACTTCATGCCAGACAACTCTGTAACGTTAGAAAGTGAAGGAAGTGGCAGCC
CAGTTTACATTACTTCATTTGGACAGGGAACACAATCTTTCAGCCTGGACCTGTGTAGAAGAAATTTGG
AGAAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001277171

Insert Size:	1458 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001277171.1</u>
RefSeq Size:	3699 bp
RefSeq ORF:	1458 bp
Locus ID:	166863
UniProt ID:	<u>Q8TBY0</u>
Cytogenetics:	4q32.1
MW:	55.3 kDa