

## Product datasheet for **SC336405**

### SLC38A1 (NM\_001278387) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** SLC38A1 (NM\_001278387) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** SLC38A1  
**Synonyms:** ATA1; NAT2; SAT1; SNAT1  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC336405 representing NM\_001278387.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

ATGATGCATTTCAAAGTGGACTCGAATTAAGTGGAGTTGCAAACATGACAGTGCCCGAGGATGATAAC
ATTAGCAATGACTCCAATGATTTACCCGAAGTAGAAAAAGGTCAGATAAATAGCAAGTTTATTTCTGAT
CGTGAAAGTAGAAGAAGTCTCACAAACAGCCATTTGGAAAAAAGAAGTGTGATGAGTATATTCCAGGT
ACAACCTCCTTAGGCATGTCTGTTTTAACCTAAGCAACGCCATTATGGGCAGTGGGATTTGGGACTC
GCCTTTGCCCTGGCAAACACTGGAATCCTACTTTTTCTGGTACTTTTGACTTCAGTGACATTGCTGTCT
ATATATCAATAAACCTCCTATTGATCTGTTCAAAGAAACAGGCTGCATGGTGTATGAAAAGCTGGGG
GAACAAGTCTTTGGCACCACAGGAAGTTCGTAATCTTTGGAGCCACCTCTCTACAGAACACTGGAGCA
ATGCTGAGCTACCTCTTCATCGTAAAAAATGAACTACCCTCTGCCATAAAGTTTCTAATGGGAAAGGAA
GAGACATTTTCAGCCTGGTACGTGGATGGCCGCTTCTGGTGGTATAGTTACCTTTGGCATAATTTCTC
CCTCTGTGTCTCTGAAGAAGTATAGGATCTTTGGCTATACTAGTGGATTTTCTTGAGCTGTATGGTT
TTTTTCTAATTGTGGTTATTTACAAGAAATTTCAAATCCCTGCATTGTTCCAGAGCTAAATTTCAACA
ATAAGTCTAATTTCAACAAATGCTGACACGTGTACGCCAAAATATGTTACCTTCAATTTCAAAGACCGTG
TATGCTTTACCCACCATTGCATTTGCATTTGTTTGCCACCCGTCAGTCCGCCAATTTACAGTGAGCTT
AAAGACCGATCACAGAAAAAATGCAGATGGTTTCAAACATCTCCTTTTTGCCATGTTTGTATGTAC
TTCTTGACTGCCATTTTGGCTACTTGACATTCTATGACAACGTGCAGTCCGACCTCCTTCACAAATAT
CAGAGTAAAGATGACATTCTCATCCTGACAGTGGCGCTGGCTGTCTGTTGCTGTGATCCTCACAGTG
CCGGTGTATTTTTACCGTTTCGTTTCATCTTTATTTGAACTGGCTAAGAAAACAAAGTTTAAATTTATGT
CGTCATACCGTGGTTACCTGCATACTCTGGTTGTTATCAACTTGTGGTGTATCATAACCTCCATG
AAGGATATTTTTGGAGTCGTAGGAGTTACATCTGCTAACATGCTTATTTTATTCTTCTTCATCTCTT
TATTTAAAAATCACAGACCAGGATGGAGATAAAGGAACTCAAAGAATTTGGGCTGCCCTTTCTTGGGC
CTGGGGTGTGTTCTCCTTGGTCAGCATTCCCTTGGTCATCTATGACTGGCCCTGCTCATCGAGTAGT
GACGAAGGCCACTGA
  
```

**Restriction Sites:** Sgfl-Mlul



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001278387

**Insert Size:** 1464 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001278387.1](#)

**RefSeq Size:** 8623 bp

**RefSeq ORF:** 1464 bp

**Locus ID:** 81539

**UniProt ID:** [Q9H2H9](#)

**Cytogenetics:** 12q13.11

**Protein Families:** Transmembrane

**MW:** 54 kDa

**Gene Summary:** Amino acid transporters play essential roles in the uptake of nutrients, production of energy, chemical metabolism, detoxification, and neurotransmitter cycling. SLC38A1 is an important transporter of glutamine, an intermediate in the detoxification of ammonia and the production of urea. Glutamine serves as a precursor for the synaptic transmitter, glutamate (Gu et al., 2001 [PubMed 11325958]).[supplied by OMIM, Mar 2008]  
Transcript Variant: This variant (3) is the longest transcript and encodes the predominant isoform (1). Variants 1, 2, 3, 4 and 5 encode the same isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.