

Product datasheet for **SC336415**

beta IV Tubulin (TUBB4A) (NM_001289127) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	beta IV Tubulin (TUBB4A) (NM_001289127) Human Untagged Clone
Tag:	Tag Free
Symbol:	TUBB4A
Synonyms:	beta-5; DYT4; TUBB4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC336415 representing NM_001289127.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGAAGGGGGGCTGCGGACCGAGAAACTGAGGGACCGTGTCCGCCGTCTCCGCCGATCTTCCACC
CTCGCCGCCCGCAGCTCCCCCGCTCGTGCCACCGCCCGCGTCCACCCTCAGCGCCACCGCCATG
CGGGAGATCGTGCACTGCAGGCCGCGCCAGTGCGGCAACCAGATCGGGGCCAAGTTTTGGAGGTTATC
AGTGACGAACATGGCATCGACCCACAGGCACATACCATGGGGACAGTGACCTGCAACTGGAGAGGATC
AACGTGTACTACAACGAGGCCACAGGAGAAATTATGTCCCAGAGCGGTGCTGGTGGACCTGGAACCC
GGCACCATGGACTCTGTCCGTTCTGGCCCTTCGGTCAGATCTTTCGGCCGGACAACCTCGTGTGGC
CAATCCGAGCCGGCAACAACCTGGGCAAAGGGGCACTACACGGAGGGCGCAGAGCTGGTGGACGCTGTC
CTGGACGTAGTCCGAAGGAGGCCGAGAGCTGCGACTGCCTTCAGGGCTCCAGCTGACCCACTCGCTG
GGGGTGGCACGGGTCCGGAATGGGCACGCTGCTCATCAGTAAGATCCGCGAGGAGTCCAGACCCGC
ATCATGAACACCTTCAGCGTGGTCCCTCGCCAAAGTGTGAGACACGGTGGTGGAGCCCTACAACGCC
ACGCTGTCTGTGACCAGCTGGTGGAGAATACGGATGAGACCTACTGCATCGACAACGAGGCACTTAC
GACATCTGTTCCGCACCTCAAGCTGACCACCCACCTACGGGGACCTCAACCACCTGGTGTCCGGCC
ACCATGAGCGGGTACCACCTGCCTGCCTTCCCGGGCCAGCTGAACGCCGACCTGCGCAAGCTGGCC
GTCAACATGGTTCCTTCTCGCTGCACTTCTCATGCCCGGCTTCGCACCCCTGACCAGCCGGGGC
AGCCAGCAGTACCGGGCCCTGACGGTGGCCGAGCTCACCCAGCAGATGTTTCGATGCCAAGAATGATG
GCGGGTGGCACCAGCCGCGCCACGGCCGCTACCTGACCGTGGCCGCGTGTCCGGGGCCGATGTCATG
AAGGAGGTGGACGAGCAGATGCTGAGCGTGCAGAGCAAGAACAGCAGCTACTTCGTGGAGTGGATCCCC
AACACGTGAAGACGGCCGTGTGCGACATCCCGCCCGCGGCTGAAGATGGCCGCGACCTCATCGGC
AACAGCACGGCCATCCAGGAGCTGTTCAAGCGCATCTCCGAGCAGTTCACGGCCATGTTCCGGCGCAAG
GCCTTCTTGCACCTGGTACACGGCGAGGCATGGACGAGATGGAGTTCACCGAGGCCGAGAGCAACATG
AATGACCTGGTATCTGAGTACCAGCAGTACCAGGACGCCACGGCCGAGGAGGGCGAGTTCGAGGAGGAG
GCGGAGGAGGAGGTGGCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001289127

Insert Size:	1470 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001289127.1
RefSeq Size:	2483 bp
RefSeq ORF:	1470 bp
Locus ID:	10382
UniProt ID:	P04350
Cytogenetics:	19p13.3
Protein Families:	Druggable Genome
Protein Pathways:	Gap junction, Pathogenic Escherichia coli infection
MW:	53.8 kDa
Gene Summary:	<p>This gene encodes a member of the beta tubulin family. Beta tubulins are one of two core protein families (alpha and beta tubulins) that heterodimerize and assemble to form microtubules. Mutations in this gene cause hypomyelinating leukodystrophy-6 and autosomal dominant torsion dystonia-4. Alternate splicing results in multiple transcript variants encoding different isoforms. A pseudogene of this gene is found on chromosome X. [provided by RefSeq, Jan 2014]</p> <p>Transcript Variant: This variant (2) has an alternate splice site in the 5' coding region, compared to variant 1. The resulting isoform (2) lacks an internal segment, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>